# Laboratory Equipment

in customized premium quality

MADE IN GERMANY

# Let's shake it

## Edmund Bühler GmbH

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## Edmund Bühler GmbH

Founded by Edmund Bühler in the university city of Tübingen in 1880, Edmund Bühler GmbH is still owned by the family today.

More than 140 years later, Edmund Bühler GmbH has developed to become a company which operates throughout the world. The production and sales of conventional laboratory devices is traditionally a main focus of the company. A second focal point is the development and production of research equipment, in particular for the quickly growing field of materials science.

In the laboratories of the world, the „Bühler Shaker“ is considered to be a synonym for the efficiency of the company. Competence, experience and the commitment to excellence in quality, which is steadily anchored in the company philosophy, form the basis for a unique program of technically mature products.

Do not expect less.

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## Miniature Shakers

KM CO<sub>2</sub>

KM 2 / KM 2 Akku

RM 2


TiMix 2

## Vario Shakers

TL 2

TL 10

WS 10

KL 2

## Compact Shakers

KS 15 A / B

KS 15 A / B control

TiMix 5

TiMix 5 control

## Universal Shakers

SM 30 A / B / C

SM 30 A / B / C control

SM 30 AT control

Beer Bottle Shaker

## Incubator Hoods

TH 15

TH 30

## Multi-Flask Shakers


VKS 75 A / B control

Barrel Shaker VKS 75 F control

Barrel Shaker VKS 75 F control

Champagne Bottle Shaker

## Accessories Loading Capacities

## Infrared Heaters

IRB 1

IRB 2

IRB 6

## Customized Solutions

A large, artistic water splash graphic on the left side of the page. It features a blue splash at the top left and a larger, more complex green splash below it, with many small droplets and bubbles trailing off to the right. The background is white with a light gray vertical band on the right.

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## Description

Stainless Steel Shaker especially for incubators with CO<sub>2</sub> concentration

## Basic equipment

Basic device incl. basic platform with rubber mat and external control unit

## Technical Data

|                                                  |                              |
|--------------------------------------------------|------------------------------|
| Motion:                                          | orbital                      |
| Shaking platform:                                | 200 x 295 mm                 |
| Max. load:                                       | 2 kg                         |
| Shaking speed:                                   | 5 - 220 rpm, in steps of 5   |
| Stroke:                                          | 20 mm                        |
| Runtime:                                         | 0 - 99 h 59 min / continuous |
| CO <sub>2</sub> concentration:                   | 20 %                         |
| Electrical supply:                               | 100 - 240 V, 50 / 60 Hz      |
| Enclosure protection:                            | IP 21                        |
| Power consumption:                               | max. 6 W                     |
| Heat emission:                                   | approx. 0.5 - 1.5 W          |
| Ambient temperature:                             | 4°C to 45°C                  |
| Relative humidity:                               | ~ 95 % (non-condensing)      |
| Length of standard cable (KM CO <sub>2</sub> ):  | 200 cm                       |
| Length of ribbon cable (KM CO <sub>2</sub> -FL): | 120 cm                       |
| Dimensions (w x d x h):                          | 230 x 340 x 145 mm           |
| Weight:                                          | 7.6 kg                       |

## Loading Capacity (Examples)

|                                              |          |
|----------------------------------------------|----------|
| Erlenmeyer flasks 100 ml                     | 9 pieces |
| Erlenmeyer flasks 500 ml                     | 2 pieces |
| Further loading possibilities see page 47/48 |          |

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The external control unit with LED display enables easy operation outside the incubator. Cable feed-through required in the incubator.

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The external control unit with LED display and ribbon cable enables easy operation outside the incubator (KM CO<sub>2</sub>-FL).

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**Universal tray KM Mini**  
No. 0052 119

Coated tray with drillings for individual loading with spring clamps or test tube racks.  
Shaking platform 200 x 295 mm.

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**Universal tray KM**  
No. 0052 030

Coated tray with drillings for individual loading with spring clamps or test tube racks.  
Shaking platform 356 x 300 mm.

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### Description

Highly reliable shaker which requires little space

### Basic equipment

Basic device incl. basic platform and rubber mat

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### Technical Data

|                         |                              |
|-------------------------|------------------------------|
| Motion:                 | orbital                      |
| Shaking platform:       | 200 x 295 mm                 |
| Max. load:              | 2 kg                         |
| Shaking speed:          | 5 - 220 rpm, in steps of 5   |
| Stroke:                 | 20 mm                        |
| Runtime:                | 0 - 99 h 59 min / continuous |
| Electrical supply:      | 100 - 240 V, 50 / 60 Hz      |
| Enclosure protection:   | IP 21                        |
| Power consumption:      | max. 6 W                     |
| Heat emission:          | approx. 0.5 - 2 W            |
| Ambient temperature:    | 4°C to 45°C                  |
| Relative humidity:      | ~ 95 % (non-condensing)      |
| Dimensions (w x d x h): | 250 x 415 x 145 mm           |
| Weight:                 | 8.1 kg                       |

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### Description

Well suited for operation in refrigerators or incubators

### Basic equipment

Basic device incl. basic platform and rubber mat, mains / charger unit

### Technical Data

|                                    |                              |
|------------------------------------|------------------------------|
| Motion:                            | orbital                      |
| Shaking platform:                  | 200 x 295 mm                 |
| Max. load:                         | 2 kg                         |
| Shaking speed:                     | 5 - 220 rpm, in steps of 5   |
| Stroke:                            | 20 mm                        |
| Charging time (battery):           | 10 hours                     |
| Runtime of battery:                | min. 40 hours                |
| Runtime with mains cable:          | 0 - 99 h 59 min / continuous |
| Electrical supply battery charger: | 100 - 240 V, 50 / 60 Hz      |
| Internal electrical supply:        | 1 battery 12 V / 7Ah         |
| Enclosure protection:              | IP 21                        |
| Power consumption:                 | max. 12 W                    |
| Heat emission:                     | approx. 0.5 - 2 W            |
| Ambient temperature:               | 4°C to 45°C                  |
| Relative humidity:                 | ~ 95 % (non-condensing)      |
| Dimensions (w x d x h):            | 250 x 415 x 145 mm           |
| Weight:                            | 10.5 kg                      |

### Loading Capacity (Examples)

|                         |           |
|-------------------------|-----------|
| Erlenmeyer flasks 10 ml | 20 pieces |
| Erlenmeyer flasks 50 ml | 12 pieces |

Further loading possibilities see page 47/48

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
### Description

Space saving tube roller with gentle tilting and rotating motion

### Basic equipment

Basic device with 6 removable plastic rolls and spill tray

### Technical Data

|                         |                              |
|-------------------------|------------------------------|
| Motion:                 | rotating, tilting motion     |
| Shaking platform:       | 200 x 300 mm                 |
| Dimensions of rolls:    | 300 x Ø 32 mm                |
| Number of rolls:        | 6 pc                         |
| Max. load:              | 2 kg                         |
| Shaking speed:          | 5 - 80 rpm, in steps of 5    |
| Angle of tilt:          | 3° fixed                     |
| Stroke at end of rolls: | 15 mm                        |
| Runtime:                | 0 - 99 h 59 min / continuous |
| Electrical supply:      | 100 - 240 V, 50 / 60 Hz      |
| Enclosure protection:   | IP 21                        |
| Power consumption:      | max. 6 W                     |
| Heat emission:          | approx. 0.5 - 2 W            |
| Ambient temperature:    | 5°C to 45°C                  |
| Relative humidity:      | ~ 85 % (non-condensing)      |
| Dimensions (w x d x h): | 244 x 465 x 115 mm           |
| Weight:                 | 7.3 kg                       |

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### Loading Capacity (Examples)

Suitable for all standard tube sizes, flasks and universal vessels up to 80 mm in diameter thanks to removable rolls

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### Overhead Rack


No. 0060 003

with 12 drillings à Ø 29,5 mm  
(other drillings on request)

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For easy cleaning the plastic rolls (no. 0060 002) as well as the spill tray (no. 0016 416) can be removed.

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### Description

Special shaker for small samples in microplates

### Basic equipment

Basic device incl. basic platform with holders for 4 microplates

### Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | exact orbital motion                                           |
| Shaking platform:       | 200 x 295 mm for 4 microplates                                 |
| Max. load:              | 2 kg                                                           |
| Shaking speed:          | 100 - 1400 rpm, steplessly adjustable                          |
| Stroke:                 | 3 mm                                                           |
| Timer:                  | 0 - 120 min / continuous                                       |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 23 W                                                      |
| Heat emission:          | approx. 6 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 220 x 410 x 110 mm                                             |
| Weight:                 | 8.5 kg                                                         |

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### Loading Capacity (Examples)

|               |                                                          |
|---------------|----------------------------------------------------------|
| 4 microplates | deepwell-plates, or other plates<br>in microplate format |
|---------------|----------------------------------------------------------|

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### Shaking platform KM No. 0052 072

For shaker TiMix 2 (instead of standard platform). Flat shaking platform with rim (20 mm) and rubber mat. Shaking platform 200 x 295 mm.

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Slow, three-dimensional motion

## Basic equipment

Basic device incl. basic tray and rubber mat

## Technical Data

|                         |                                   |
|-------------------------|-----------------------------------|
| Motion:                 | three-dimensional orbital motion  |
| Shaking platform:       | 220 x 220 mm                      |
| Max. load:              | 2 kg                              |
| Shaking speed:          | 1 - 22 rpm, steplessly adjustable |
| Angle of tilt:          | 8° fixed                          |
| Runtime:                | 0 - 99 h 59 min / continuous      |
| Electrical supply:      | 100 - 240 V, 50 / 60 Hz           |
| Enclosure protection:   | IP 21                             |
| Power consumption:      | max. 6 W                          |
| Heat emission:          | approx. 0.5 - 2 W                 |
| Ambient temperature:    | 5°C to 45°C                       |
| Relative humidity:      | ~ 85 % (non-condensing)           |
| Dimensions (w x d x h): | 250 x 415 x 118 mm                |
| Weight:                 | 9.1 kg                            |

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
Slow, three-dimensional motion

## Basic equipment

Basic device incl. basic tray and rubber mat

## Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | three-dimensional orbital motion                               |
| Shaking platform:       | 340 x 340 mm                                                   |
| Max. load:              | 10 kg                                                          |
| Shaking speed:          | 1 - 23 rpm, steplessly adjustable                              |
| Angle of tilt:          | 0° to 9° adjustable                                            |
| Timer:                  | 0 - 120 min / continuous                                       |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 18 W                                                      |
| Heat emission:          | approx. 6 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 355 x 450 x 205 mm                                             |
| Weight:                 | 18.5 kg                                                        |

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### Multi-storey rack system

No. 0051 247

With 2 trays (340 x 340 mm) and rubber mat.  
For shakers WS 10 and TL 10.  
Height: 200 mm

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(340 x 340 mm)

No. 0051 502

Additional tray for multi-storey rack WS / TL incl. rubber mat.

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Slow, linear tilting motion

## Basic equipment

Basic device incl. basic tray and rubber mat

## Technical Data

|                         |                                   |
|-------------------------|-----------------------------------|
| Motion:                 | linear tilting motion             |
| Shaking platform:       | 340 x 340 mm                      |
| Max. load:              | 10 kg                             |
| Shaking speed:          | 1 - 23 rpm, steplessly adjustable |
| Angle of tilt:          | 5° fixed                          |
| Timer:                  | 0 - 120 min / continuous          |
| Electrical supply:      | 115 - 230 V, 50 / 60 Hz           |
| Enclosure protection:   | IP 21                             |
| Power consumption:      | max. 8 W                          |
| Heat emission:          | approx. 3 W                       |
| Ambient temperature:    | 5°C to 50°C                       |
| Relative humidity:      | ~ 85 %                            |
| Dimensions (w x d x h): | 355 x 450 x 205 mm                |
| Weight:                 | 18.5 kg                           |

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The „2 in 1“ shaker

## Basic equipment

Basic device incl. rack system with 3 clamping strips h and rubber mat

## Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | can be switched from orbital to reciprocating (horizontal)     |
| Shaking platform:       | 300 x 300 mm                                                   |
| Max. load:              | 2 kg                                                           |
| Shaking speed:          | 30 - 420 rpm, steplessly adjustable                            |
| Stroke:                 | 8 mm                                                           |
| Timer:                  | 0 - 120 min / continuous                                       |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 23 W                                                      |
| Heat emission:          | approx. 4 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 355 x 455 x 195 mm                                             |
| Weight:                 | 13.5 kg                                                        |

## Loading Capacity (Examples)

|                                              |          |
|----------------------------------------------|----------|
| Erlenmeyer flasks 250 ml                     | 9 pieces |
| Separating funnels 250 ml                    | 1 pieces |
| Further loading possibilities see page 47/48 |          |

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Nr. / No. 6173 000

## Description

Reliable, robust shaker for daily use

## Basic equipment

Basic device without shaking plate, without rack system

## Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion KS 15 A:         | orbital                                                        |
| Motion KS 15 B:         | reciprocating (horizontal)                                     |
| Shaking platform:       | 400 x 300 mm                                                   |
| Max. load:              | 15 kg                                                          |
| Shaking speed:          | 30 - 420 rpm, steplessly adjustable                            |
| Stroke:                 | 17 mm                                                          |
| Timer:                  | 0 - 120 min / continuous                                       |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 25 W                                                      |
| Heat emission:          | approx. 5 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 510 x 490 x 150 mm                                             |
| Weight:                 | 18 kg                                                          |

## Loading Capacity (Examples)

|                                              |           |
|----------------------------------------------|-----------|
| Erlenmeyer flasks 100 ml                     | 20 pieces |
| Erlenmeyer flasks 250 ml                     | 12 pieces |
| Erlenmeyer flasks 500 ml                     | 6 pieces  |
| Erlenmeyer flasks 1000 ml                    | 4 pieces  |
| Separating funnels 250 ml                    | 1 pieces  |
| Further loading possibilities see page 47/48 |           |

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### Rack system Combifix KS

No. 0052 071

Rack system with 3 clamping strips h and rubber mat.  
Shaking platform 400 x 275 mm, height approx. 80 mm.  
For fastening upright vessels.  
With additional strips also suitable for other vessels.

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### Universal tray KS

No. 0051 471

Coated tray with drillings for individual fastening of spring clamps or test tube racks.  
Shaking platform 430 x 300 mm.

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### Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion KS 15 A control: | orbital                                                        |
| Motion KS 15 B control: | reciprocating (horizontal)                                     |
| Shaking platform:       | 400 x 300 mm                                                   |
| Max. load:              | 15 kg                                                          |
| Shaking speed:          | 30 - 420 rpm, in steps of 5                                    |
| Stroke:                 | 17 mm                                                          |
| Runtime:                | programmable / continuous                                      |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 25 W                                                      |
| Heat emission:          | approx. 7 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 510 x 490 x 150 mm                                             |
| Weight:                 | 19 kg                                                          |

### Loading Capacity (Examples)

|                                              |           |
|----------------------------------------------|-----------|
| Erlenmeyer flasks 100 ml                     | 20 pieces |
| Erlenmeyer flasks 250 ml                     | 12 pieces |
| Erlenmeyer flasks 500 ml                     | 6 pieces  |
| Erlenmeyer flasks 1000 ml                    | 4 pieces  |
| Separating funnels 250 ml                    | 1 pieces  |
| Further loading possibilities see page 47/48 |           |

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#### Defined stop\*

Nr. 0060 019

Stop always in the same position (KS 15 A control; KS 15 B control on request)

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### Description

High speed for optimal mixing of small sample volumes

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
### Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | exact orbital motion                                           |
| Shaking platform:       | 400 x 300 mm for 8 microplates                                 |
| Max. load:              | 5 kg                                                           |
| Shaking speed:          | 100 - 1400 rpm, steplessly adjustable                          |
| Stroke:                 | 3 mm                                                           |
| Timer:                  | 0 - 120 min / continuous                                       |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 30 W                                                      |
| Heat emission:          | approx. 9 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 510 x 490 x 150 mm                                             |
| Weight:                 | 19 kg                                                          |

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**Standard rack system**  
No. 0052 101

For 8 standard microplates, microwell or deepwell plates, or other plates of the same size.

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**Rack system Combifix KS**  
No. 0052 071

Rack system with 3 clamping strips h and rubber mat. Shaking platform 400 x 275 mm, height approx. 80 mm. For fastening upright vessels. With additional strips also suitable for other vessels.

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**Additional tray**  
No. 0052 102

As a second storey on the standard rack system TiMix 5. Distance between the 2 trays max. 55 mm.

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**Universal tray KS**  
No. 0051 471

Coated tray with drillings for individual fastening of spring clamps or test tube racks. Shaking platform 430 x 300 mm.

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**Rack system with high clamping pins**  
No. 0052 104

For max. 24 standard microplates, placed directly on top of each other.

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
**Rack system with metal clamping pins**  
No. 0052 125

For max. 48 standard microplates, placed directly on top of each other.

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### Description

Top technology - also for integration in analytical robot systems

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### Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | exact orbital motion                                           |
| Shaking platform:       | 400 x 300 mm for 8 microplates                                 |
| Max. load:              | 5 kg                                                           |
| Shaking speed:          | 100 - 1400 rpm, in steps of 10                                 |
| Stroke:                 | 3 mm<br>(Special stroke on request)                            |
| Runtime:                | programmable / continuous                                      |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 30 W                                                      |
| Heat emission:          | approx. 11 W                                                   |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 510 x 490 x 150 mm                                             |
| Weight:                 | 20 kg                                                          |

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Nr. / No. 6105 000

## Description

Table-top shaker for high loads and continuous operation

## Basic equipment

Basic device with shaking plate and rubber mat, without rack system

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**Rack system  
Combifix SM A**  
No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels.

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**2-storey top frame SM**  
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 332 mm. Distance between the 2 storeys: 300 mm (height sufficient for 1000 ml Erlenmeyer flasks).

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**Rack system  
Combifix SM B**  
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v. For fastening horizontal vessels between the clamping strips, or as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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**Rack system  
Combifix SM C**  
No. 0050 156

Rack system with 2 clamping strips h, 2 clamp strips and 1 spring strip. For fastening separating funnels. Standard version of clamp strip: 5 clamps NS 19 (D20) and 2 clamps NS 29 (D32). With additional strips also suitable for other vessels.

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**Universal tray SM**  
No. 0051 472

Coated tray with drillings for individual fastening of spring clamps or test tube racks. Shaking platform 560 x 400 mm.

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## Sliding Plates

Sliding plates with Combifix SM (rack systems Combifix SM, but with additional base plate).

**Sliding Plate  
with Combifix SM A**  
No. 0051 484

**with Combifix SM B**  
No. 0051 485

**with Combifix SM C**  
No. 0051 486

**Mounting set**  
No. 0052 074

for universal trays, or sliding plates on the 2-storey top frame SM.

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## Description

Table-top shaker with programmable shaking tasks for the mixing of difficult media

## Basic equipment

Basic device with shaking plate and rubber mat, without rack system.  
Programmable control with touch display  
(Technical details of control see page 44)

## Technical Data

|                         |                                                                       |
|-------------------------|-----------------------------------------------------------------------|
| Motion SM 30 A control: | orbital                                                               |
| Motion SM 30 B control: | reciprocating (horizontal)                                            |
| Motion SM 30 C control: | can be switched<br>from orbital to reciprocating                      |
| Shaking platform:       | 560 x 400 mm                                                          |
| Max. load:              | 30 kg                                                                 |
| Shaking speed:          | 15 - 300 rpm, in steps of 5                                           |
| Stroke SM 30 A control: | 26 mm                                                                 |
| Stroke SM 30 B control: | 30 mm<br>46 mm (Option: No. 0052 212)<br>50 mm (Option: No. 0052 216) |
| Stroke SM 30 C control: | 26 mm                                                                 |
| Runtime:                | programmable / continuous                                             |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order        |
| Enclosure protection:   | IP 21                                                                 |
| Power consumption:      | max. 46 W                                                             |
| Heat emission:          | approx. 8 W                                                           |
| Ambient temperature:    | 5°C to 50°C                                                           |
| Relative humidity:      | ~ 85 %                                                                |
| Dimensions (w x d x h): | 680 x 610 x 160 mm                                                    |
| Weight:                 | 33 kg                                                                 |

## Loading Capacity (Examples)

|                                              |           |
|----------------------------------------------|-----------|
| Erlenmeyer flasks 100 ml                     | 42 pieces |
| Erlenmeyer flasks 250 / 300 ml               | 20 pieces |
| Separating funnels 250 ml                    | 6 pieces  |
| Separating funnels 1000 ml                   | 2 pieces  |
| Test tube racks                              | 4 pieces  |
| Further loading possibilities see page 47/48 |           |

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### Defined stop\*

No. 0060 020

Stop always in the same position (SM 30 A control and SM 30 B control)

\*Installation only possible in our works

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**2-storey top frame SM**  
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 332 mm. Distance between the 2 storeys: 300 mm. (height sufficient for 1000 ml Erlenmeyer flasks)

Alternatively, the rack systems Combifix SM or the universal tray SM can be mounted on the 2-storey top frame SM.

For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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### Description

Table-top shaker with heating plate and programmable shaking tasks for mixing of difficult media

### Basic equipment

Basic device with heating plate, without rack system.  
Programmable control with touch display  
(Technical details of control see page 44)

### Technical Data

|                         |                                               |
|-------------------------|-----------------------------------------------|
| Motion:                 | orbital                                       |
| Shaking platform:       | 560 x 400 mm                                  |
| Max. load:              | 25 kg                                         |
| Shaking speed:          | 15 - 300 rpm, in steps of 5                   |
| Stroke:                 | 26 mm                                         |
| Runtime:                | programmable / continuous                     |
| Temperature range:      | up to 120°C                                   |
| Control accuracy:       | approx. +/- 1°C                               |
| Electrical supply:      | 230 V, 50 / 60 Hz<br>115 V version on request |
| Enclosure protection:   | IP 21                                         |
| Power consumption:      | max. 450 W                                    |
| Heat emission:          | 400 W                                         |
| Ambient temperature:    | 5°C to 50°C                                   |
| Relative humidity:      | ~ 85 %                                        |
| Dimensions (w x d x h): | 680 x 610 x 160 mm                            |
| Weight:                 | 40 kg                                         |

### Loading Capacity (Examples)

|                                              |           |
|----------------------------------------------|-----------|
| Erlenmeyer flasks 100 ml                     | 42 pieces |
| Erlenmeyer flasks 250 / 300 ml               | 20 pieces |
| Separating funnels 250 ml                    | 6 pieces  |
| Separating funnels 1000 ml                   | 2 pieces  |
| Further loading possibilities see page 47/48 |           |

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**Rack system  
Combifix SM A**  
No. 0050 154

Basic rack with 5 clamping strips h, 2 clamp strips and 1 spring strip. For fastening different upright vessels. With additional strips also suitable for other vessels.

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**Rack system  
Combifix SM B**  
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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The below specified shaking speeds against load are approximate values. Depending on specific properties of substances to be shaken these values can differ marginally.

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#### Description

Programmable shaking tasks for the mixing of difficult media

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#### Technical Data

|                         |                                                                |
|-------------------------|----------------------------------------------------------------|
| Motion:                 | reciprocating (horizontal)                                     |
| Shaking platform:       | 560 x 400 mm                                                   |
| Max. load:              | 30 kg                                                          |
| Shaking speed:          | 15 - 240 rpm, in steps of 5                                    |
| Stroke:                 | 50 mm for Gushing Test<br>(No. 0052 216)                       |
| Runtime:                | programmable / continuous                                      |
| Electrical supply:      | 230 V or 115 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                          |
| Power consumption:      | max. 46 W                                                      |
| Heat emission:          | approx. 8 W                                                    |
| Ambient temperature:    | 5°C to 50°C                                                    |
| Relative humidity:      | ~ 85 %                                                         |
| Dimensions (w x d x h): | 680 x 610 x 160 mm                                             |
| Weight:                 | 33 kg                                                          |

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2 x Spannleiste h Nr. 0050 400  
1 x Spannleiste q Nr. 0050 403

## 1-storey rack system No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels. (see picture)

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## 2-storey top frame No. 0052 065

Top frame with 2 coated trays with drillings for mounting rack systems Combifix SM or for spring clamps or test tube racks.

Total height: 332 mm;  
Distance between the 2 storeys: 300 mm.

Dimensions (w x d x h):  
680 x 420 x 480 mm  
Weight: 7.7 kg

Alternatively, the universal tray SM can be mounted on the 2-storey top frame SM. For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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## 3-storey top frame No. 0039 206

Top frame with 3 coated trays for mounting the rack systems Combifix SM. Total height: 400 mm; Distance between the storeys: 177 mm.

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### Maximum Shaking Speed Beer Bottle Shaker SM 30 B

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The below specified shaking speeds against load are approximate values.

Depending on specific properties of substances to be shaken these values can differ marginally.

| Typ / Type | Aufsatz / Rack system | Hub / stroke (mm) | 75 rpm | 125 rpm | 180 rpm | 200 rpm | 240 rpm |
|------------|-----------------------|-------------------|--------|---------|---------|---------|---------|
| SM 30 B    | 1-etagig / 1-storey   | 50 mm             | 30 kg  | 30 kg   | 30 kg   | 20 kg   | 5 kg    |
| SM 30 B    | 2-etagig / 2-storey   | 50 mm             | 30 kg  | 20 kg   | 15 kg   | 5 kg    |         |
| SM 30 B    | 3-etagig / 3-storey   | 50 mm             | 30 kg  | 20 kg   |         |         |         |

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Nr. / No. 6161 000

**Description**

Constant temperatures; also for high vessels  
(e.g. 1000 ml Erlenmeyer flasks)  
Can be combined with all compact shakers or used as independent unit

**Basic equipment**

Coated metal housing with 3 windows, the front door opens upwards

**Technical Data**

Temperature range: 5°C above ambient to +50°C  
Temperature accuracy: +/- 1% of the set value  
Inner dimensions (w x d x h): 490 x 410 x 310 mm  
(height sufficient for 1000 ml Erlenmeyer flasks)  
Temperature sensor: PT 100  
Air circulation: 240 m³/h  
Electrical supply: 230 V, 50 / 60 Hz  
115 V version on request  
Power consumption: max. 345 W  
External Dimensions (w x d x h): 510 x 450 x 455 mm  
Weight: 21 kg  
Safety cut-off of the shaker if combined with a Bühler Compact Shaker  
Standard temperature controller 5°C above ambient to +50°C

**Option**

**Extended temperature range up to +60°C**

No. 0052 105

Installation of an additional heater in order to reach temperatures up to max. +60°C. Installation only possible in our works.

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**Shelf TH 15**

No. 0052 126

Additional tray for samples which  
have to be temperature controlled  
but need no shaking.  
Dimensions: 378 x 318 mm.

TH 30 mit SM 30 control, Universaltray und Zubehör  
TH 30 with SM 30 control, universal tray and accessories

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Nr. / No. 6162 000

## Description

Excellent system thanks to reproducible ambient temperatures and useful accessories  
Can be combined with all universal shakers or used as independent unit

## Basic equipment

Coated metal housing with 3 windows, the front door opens upwards, socket inside

## Technical Data

|                                  |                                                                         |
|----------------------------------|-------------------------------------------------------------------------|
| Temperature range:               | 5°C above ambient to +50°C                                              |
| Temperature accuracy:            | +/- 1% of the set value                                                 |
| Inner dimensions (w x d x h):    | 660 x 540 x 430 mm<br>(height sufficient for 2000 ml Erlenmeyer flasks) |
| Lighting:                        | standard light 15 W / UV light 15 W                                     |
| Temperature sensor:              | PT 100                                                                  |
| Air circulation:                 | 240 m³/h                                                                |
| Electrical supply:               | 230 V, 50 / 60 Hz<br>115 V version on request                           |
| Power consumption:               | max. 620 W                                                              |
| External Dimensions (w x d x h): | 680 x 610 x 560 mm                                                      |
| Weight:                          | 33 kg                                                                   |

Safety cut-off of UV light and stop of the shaker when the front door is opened if combined with a Bühler Universal Shaker  
Standard temperature controller 5°C above ambient to +50°C

## Options

### Extended temperature range up to +60°C

No. 0052 105

Installation of an additional heater in order to reach temperatures up to max. +60°C. Installation only possible in our works.

### Cooling device TH 30

No. 0051 581

Installation of a cooling coil and connectors for an external cooling system (refrigerated circulator or flow-through cooler). Depending on the cooling capacity of the external cooler temperatures of approx. 5°C below room temperature can be reached (300 - 500 W required). Installation only possible in our works.

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**Shelf TH 30**  
No. 0052 123

Additional tray for samples which have to be temperature controlled but need no shaking.  
Dimensions: 365 x 558 mm.

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**2-storey top frame SM / TH**  
No. 0052 117

2-storey top frame for all SM shakers. Dimensions adapted to TH 30. For vessels with a max. height of 14 cm (Erlenmeyer flasks 200 / 250 ml).  
Dimensions (w x d x h): 600 x 402 x 192 mm.  
For fastening the vessels, either universal trays SM or sliding plates with Combifix SM can be used.

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## Technical Data

Motion VKS 75 A control: orbital  
Motion VKS 75 B control: reciprocating (horizontal)  
Shaking platform: 760 x 600 mm  
Max. load: 75 kg  
Shaking speed: 30 - 200 rpm, in steps of 5  
Stroke VKS 75 A control: 26 mm  
Stroke VKS 75 B control: 50 mm  
36 mm (Option: No. 0052 213)  
60 mm (Option: No. 0052 214)  
80 mm (Option: No. 0052 215)  
programmable / continuous  
Runtime:  
Electrical supply: 100 - 240 V, 50 - 60 Hz  
please indicate in case of order  
Enclosure protection: IP 21  
Power consumption: max. 70 W  
Heat emission: approx. 15 W  
Ambient temperature: 5°C to 50°C  
Relative humidity: ~ 85 %  
Dimensions (w x d x h): 1050 x 835 x 250 mm  
Weight: 110 kg

## Loading Capacity (Examples)

Erlenmeyer flasks 500 ml 32 pieces  
Erlenmeyer flasks 2000 ml 12 pieces  
Separating funnels 1000 ml 6 pieces  
Further loading possibilities see page 47/48

## Options

### USB 2.0 interface

No. 0052 931

for the communication between PC and shaker

### Control interface

No. 0052 932

for the integration of the shaker in robotic systems (SPS)

### Safety Package\*

No. 0060 025

Monitoring of motor and shaker speed

### RS 232 interface

No. 0052 933

for the communication between robotic system and shaker

### Defined stop\*

Nr. 0060 021

Stop always in the same position (VKS 75 A control; VKS 75 B control on request)

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## Shaking plate VKS (1-storey)

No. 0052 070

Base plate 760 x 600 mm incl. PE plate. Required for mounting rack systems Combifix VKS or universal tray VKS.

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## Universal tray VKS (1-storey)

No. 0051 474

Coated tray with drillings for spring clamps or test tube racks. Shaking platform 760 x 600 mm. To be mounted on the shaking plate.

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## Rack system Combifix VKS A

No. 0051 487

Basic rack system with 9 clamping strips h. For secure fastening of different upright vessels. With additional strips also suitable for other vessels. To be mounted on the shaking plate.

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## Rack system Combifix VKS B

No. 0051 488

Basic rack system with 4 clamping strips h and 4 clamping strips v. For fastening horizontal vessels between the clamping strips, or as added support for high vessels. The maximum distance between the clamping strips h and v is 60 mm. To be mounted on the shaking plate.

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## Rack system Combifix VKS C

No. 0051 489

Basic rack system with 2 clamping strips h, 2 clamp strips and 1 spring strip. For fastening separating funnels. Standard version of clamp strip: 7 clamps NS 29 (D32). With additional strips also suitable for other vessels. To be mounted on the shaking plate.

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Total height: 510 mm  
Distance  
between levels: 387 mm  
Weight: 36 kg

Height sufficient for  
5000 ml Erlenmeyer flasks

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Total height: 910 mm  
Distance  
between levels: 340 mm  
Weight: 43 kg

Height sufficient for  
3000 ml Erlenmeyer flasks

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## Platforms for 2- and 3-storey rack systems VKS

### Sliding Plate with Combifix VKS A

No. 0051 493

Base plate 760 x 600 mm with basic  
rack and 9 clamping strips h.

### with Combifix VKS B

No. 0051 494

Base plate 760 x 600 mm with basic  
rack, 4 clamping strips h and 4 clam-  
ping strips v.

### with Combifix VKS C

No. 0051 495

Base plate 760 x 600 mm with basic  
rack and 2 clamping strips h,  
2 cramp strips and 1 spring strip.

### Universal tray VKS

No. 0051 496

Coated tray with drillings for spring  
clamps or test tube racks.  
Shaking platform 760 x 600 mm.

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## Alternative platforms for 2- and 3-storey rack systems VKS

### Sliding Plate with Combifix SM A

No. 0051 484

Base plate 560 x 400 mm with basic  
rack and 5 clamping strips h.

### with Combifix SM B

No. 0051 485

Base plate 560 x 400 mm with basic  
rack, 4 clamping strips h and 4 clam-  
ping strips v.

### with Combifix SM C

No. 0051 486

Base plate 560 x 400 mm with basic  
rack and 2 clamping strips h,  
2 cramp strips and 1 spring strip.

### Universal tray SM

No. 0051 472

Coated tray with drillings for spring  
clamps or test tube racks.  
Shaking platform 560 x 400 mm.

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The below specified shaking speeds against load are approximate values.  
Depending on specific properties of substances to be shaken these values can differ marginally.

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#### Technical Data

|                         |                                                             |
|-------------------------|-------------------------------------------------------------|
| Motion:                 | reciprocating (horizontal)                                  |
| Shaking platform:       | 760 x 600 mm                                                |
| Max. load:              | 75 kg                                                       |
| Shaking speed:          | up to 80 rpm, in steps of 5                                 |
| Stroke optional:        | 120 - 150 mm in 10 mm steps                                 |
| Runtime:                | programmable / continuous                                   |
| Electrical supply:      | 100 - 240 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:   | IP 21                                                       |
| Power consumption:      | max. 90 W                                                   |
| Heat emission:          | approx. 15 W                                                |
| Ambient temperature:    | 5°C to 50°C                                                 |
| Relative humidity:      | ~ 85 %                                                      |
| Dimensions (w x d x h): | 1050 x 835 x 250 mm                                         |
| Weight:                 | 110 kg                                                      |

#### Loading Capacity

Horizontal or vertical barrels up to 80 l capacity

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**Shaking plate VKS  
(1-storey)**  
No. 0052 070

Base plate 760 x 600 mm incl. PE  
plate. Required for mounting the  
rack system no. 0039 001.

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**Special Rack System for barrels**  
No. 0039 001

Basic rack with 3 fixing elements.  
For fastening horizontal or vertical  
barrels up to 80 l capacity.  
Dimensions (w x d x h):  
780 x 565 x 450 mm  
Weight: 12 kg

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**Special Rack System  
for horizontal barrels**  
No. 0039 002

Rack system with base plate.  
For fastening horizontal barrels  
up to 50 l capacity  
Dimensions (w x d x h):  
800 x 600 x 200 mm  
Weight: 15 kg


### Maximale Schüttelfrequenzen für Fass-Schüttler VKS 75 F control Maximum Shaking speed for Barrel shaker VKS 75 F control

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#### Description

Big heavy-duty shaker for high loads also during continuous operation  
(floor installation only)

#### Basic equipment

Basic device without shaking plate, without rack system.  
Programmable control with touch display  
(Technical details of control see page 44)

#### Technical Data

|                                         |                                                             |
|-----------------------------------------|-------------------------------------------------------------|
| Motion:                                 | reciprocating (horizontal)                                  |
| Shaking platform:                       | 760 x 600 mm                                                |
| Max. load:                              | 75 kg                                                       |
| Shaking speed:                          | up to 80 rpm, in steps of 5                                 |
| Stroke:                                 | 130 mm                                                      |
| Runtime:                                | programmable / continuous                                   |
| Electrical supply:                      | 100 - 240 V, 50 / 60 Hz<br>please indicate in case of order |
| Enclosure protection:                   | IP 21                                                       |
| Power consumption:                      | max. 90 W                                                   |
| Heat emission:                          | approx. 15 W                                                |
| Ambient temperature:                    | 5°C to 50°C                                                 |
| Relative humidity:                      | ~ 85 %                                                      |
| Dimensions shaker (w x d x h):          | 1050 x 835 x 250 mm                                         |
| Weight:                                 | 110 kg                                                      |
| Total weight<br>(Shaker + rack system): | 135 kg                                                      |

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##### USB 2.0 interface

No. 0052 931

for the communication between PC  
and shaker

##### Control interface

No. 0052 932

for the integration of the shaker in  
robotic systems (SPS)

##### Safety Package\*

No. 0060 025

Monitoring of motor and shaker  
speed


##### RS 232 interface

No. 0052 933

for the communication between  
robotic system and shaker

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750 x 650 x 370 mm

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Gewicht: 25 kg

**Rack System for  
champagne bottles**  
**No.0052 287**

Rack system for fastening up to 36 horizontal champagne bottles.

To show the fixation of the champagne bottles on the rack system, side walls and lid are omitted on the illustration. The rack system is completely closed.

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- Intuitive operation via capacitive touch display for easy operation also with gloves
- Non-reflecting LCD colour display 110 x 67 mm with resolution 800 x 480 dot
- Bright LED background lighting (400cd/m<sup>2</sup>)
- Display of date and real time
- Sensor for the detection of an incubator hood on the shaker
- Manual operation or automatic operation
- Pause function for the shaking tasks
- Logfile function for the recording of interruptions of the shaking process
- Possibility to change the rotational direction for optimal mixing (orbital motion)
- Automatic operation:
  - 10 programs with 10 segments each, freely definable
  - Optical and acoustic signal at the end of a program
  - Possibility to program a defined time for starting the shaking action
  - Programmable restart after power failure
  - Alpha-numeric program names

## Options

**USB 2.0 interface (socket type B)** No. 0052 931  
**for the communication between PC and shaker**  
 Plug-and-play connection with standard terminal programs.  
 Complete remote control of the shaker.

**Interface RS 232** No. 0052 933  
**for the communication between robotic system and shaker**  
 Serial interface for industrial and robotic systems.  
 Shaker completely remote controlled.


**Control interface** No. 0052 932  
**for the integration of the shaker in robotic systems (PLC)**  
 - Relais output: shaker is running / stops  
 - Relais output: signal alarm / trouble-free operation  
 - Relais output: ready for operation  
 - Optocoupler input for external start/stop  
 - Analog input for setpoint of the speed  
 - Analog output for actual speed value

**Defined Stop\***  
 Predefined stop always in the same position. Especially suited for the use of the device in robotic systems.  
 (TiMix 5 control, KS 15 A control, SM 30 A control, SM 30 B control, VKS 75 A control; KS 15 B control and VKS 75 B control on request)

**Safety Package\***  
 The speed of both motor and shaker is continuously monitored. In case of deviations, e.g. due to a defective drive belt, a warning message is displayed and an acoustic signal is audible.

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### Clamping strip h

Holder with sponge rubber cover for fastening upright vessels.

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### Clamping strip v

Holder with sponge rubber cover in combination with clamping strip h.

For fastening horizontal or high vessels.

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### Cramp strip

For fastening horizontal separating funnels.

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### Spring strip

For securing the stoppers of separating funnels or other horizontal vessels.

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For universal trays. The sizes are related to Erlenmeyer flasks, but are also suitable for round-bottom flasks, beakers, etc.

|          |          |          |          |          |          |          |          |          |          |
|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|
| 10 ml    | 25 ml    | 50 ml    | 100 ml   | 250 ml   | 500 ml   | 1000 ml  | 2000 ml  | 3000 ml  | 5000 ml  |
| 0009 642 | 0009 643 | 0009 644 | 0009 645 | 0009 646 | 0009 647 | 0009 648 | 0009 649 | 0009 653 | 0009 652 |

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To be fastened either with clamping strips on rack systems or with a hinged foot on universal trays. Dimensions 265 x 100 mm each.

With the hinged foot it is possible to vary the angle of inclination of the test tubes.

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Abmessung 289 x 100 mm.

Required for universal trays. Size fitting for all test tube racks. Dimensions 289 x 100 mm.

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<sup>1)</sup> with 5 clamping strips h (0052 115) additional <sup>2)</sup> with 4 clamping strips h (0052 115) additional <sup>3)</sup> with 3 clamping strips h (0052 115) additional

<sup>4)</sup> with 2 clamping strips h (0052 115) additional <sup>5)</sup> with 2 clamping strips h (0050 118) additional <sup>6)</sup> with 1 clamping strip h (0050 118) additional

<sup>7)</sup> with 3 clamping strips h (0050 118) additional <sup>8)</sup> with 3 clamping strips h (0050 400) additional <sup>9)</sup> with 1 clamping strip h (0050 400) additional

<sup>10)</sup> with 5 clamping strips h (0050 387) additional <sup>11)</sup> with 1 clamping strip h (0050 387) additional

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<sup>12)</sup> with 1 cramp strip (0050 206) and 1 spring strip (0050 207) additional <sup>13)</sup> Number of pieces depends on shape and dimensions of the separating funnels

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
Nr. / No. 6060 000

## Description

Compact laboratory heater with fast heating-up rate

## Basic equipment

Basic device with support cross

## Technical Data

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Max. temperature:         | 700°C                                         |
| Heating power:            | 250 W                                         |
| Heating area:             | 60 x 60 mm                                    |
| Max. ambient temperature: | 30°C                                          |
| Electrical supply:        | 230 V, 50 / 60 Hz<br>115 V version on request |
| Dimensions (w x d x h):   | 100 x 105 x 105 mm                            |
| Weight:                   | 0.7 kg                                        |

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## Power control LR 1 No. 6070 000

Stepless control of the current  
consumption between 0 - 100 %

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## Support rod for stand connection

### No. 0012 164

Stainless steel rod for integration  
of the heater into a support stand,  
length approx. 13 cm

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#### Description

Laboratory heater with power control for larger volumes

#### Basic equipment

Basic device incl. integrated power control, support cross

#### Technical Data

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Max. temperature:         | 900°C                                         |
| Heating power:            | 800 W                                         |
| Heating area:             | 100 x 100 mm                                  |
| Power control:            | stepless control between 0 - 100 %            |
| Max. ambient temperature: | 30°C                                          |
| Electrical supply:        | 230 V, 50 / 60 Hz<br>115 V version on request |
| Dimensions (w x d x h):   | 150 x 150 x 170 mm                            |
| Weight:                   | 2.5 kg                                        |

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#### Support rod for stand connection

##### No. 0012 164

Stainless steel rod for integration of the heater into a support stand, length approx. 13 cm

IRB 6 mit Stativgestell, PVC-Schiene und Halteklammern  
IRB 6 with support base, horizontal PVC rod and clamps for PVC rod

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Nr. / No. 9066 000

## Description

Modular laboratory heater with 6 x IRB 1 in a housing with 6 power control units  
Each heater can be individually used and controlled


## Basic equipment

6 x IRB 1 with support crosses, 1 x power control unit LR 6

## Technical Data

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Max. temperature:         | each 700°C                                    |
| Heating power:            | 6 x 250 W                                     |
| Heating area:             | 6 x 60 x 60 mm                                |
| Power control:            | 6 stepless controls between 0 - 100 %         |
| Max. ambient temperature: | 30°C                                          |
| Electrical supply:        | 230 V, 50 / 60 Hz<br>115 V version on request |
| Heating power:            | max. 1200 W                                   |
| Dimensions (w x d x h):   | 697 x 305 x 122 mm                            |
| Weight:                   | 9 kg                                          |

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### Support base No. 0052 060

Base with 6 support rod (stainless steel), length of the rods approx. 75 cm (without PVC rod and clamps for vessels)

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### Horizontal PVC rod for support base No. 0052 091

For fastening clamps for glass vessels

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
## Examples of customized solutions

With our comprehensive laboratory device program we can offer coherent answers for most demands which arise from practical laboratory work. We believe that our job is even more. A special challenge are the very particular areas of use for which the standard is not enough. Here we design and produce customized solutions. Customer-oriented, quality conscious and tailor made. Some examples of our customized products:

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