EURO-SCIENTIFIC bvba
LERENVELD 12
BE-2547 LINT
TEL. +32 3 489 32 42
FAX +32 3 489 32 73
E-mail info@euro-scientific.be
www.euro-scientific.be

EURO-SCIENTIFIC

Heat
is our
Passion

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

**CERAN® Hot Plates
with thermostatic
temperature control 50..500 °C**

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Glass Ceramic Hotplate CT 10

**High power, rapid, glass ceramic
hotplate up to 500°C**

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- ✓ Rapid heating
- ✓ Reasonably priced
- ✓ Energy saving

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**High power rapid hotplate with a glass
ceramic heating surface for temper-
atures up to 500°C**

The glass ceramic used for the heating surface is manufactured by the Schott company and is in use in millions of households. This glass ceramic is highly resistant to breakage and sudden changes of temperature (a temperature shock of over 700°C can be withstood without damage), it is also distortion and pore-free as well as being permeable to ultra-violet light and resistant against chemicals.

Construction: The low-density, glass ceramic hotplate, electrically heated centrally over an area of 145 mm Ø, is mounted on a high-grade, stainless steel housing.



The variable temperature control system is regulated by an energy regulator 10..100% for temperatures up to 500°C. The energy regulator is housed in a separate, heat resistant container. Thereby ensuring that the regulator is not influenced by excess heat during non-stop operation.

A lighted control lamp shows when the hotplate is switched on.

The appliance has a integrated temperature sensor to safeguard against overheating.

The ready-to-connect hotplate is delivered complete with a 1.8 m mains cable.

Best-Nr.	CT 10	Order No.
Spannung	230 - 240 V, 50-60 Hz	voltage
Leistung	1200 Watt	performance
Glaskeramik	175 mm x 175 mm	glass ceramic
Beheizte Fläche	145 mmØ	heating surface
Breite x Tiefe	200 mm x 310 mm	width and depth
Höhe	95 mm	height
Gewicht	2,8 kg	weight

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The right to make modifications serving further technical development, is reserved

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Further Fabrications

High temperature TITANIUM hot plate up to 600°C - CERAN heating plates to 5200C - Electronic DIGITAL hot plates to 450°C - Precision hot plates - Stand hot plates with supporting rod - Multi-purpose, extremely deep, heating bath - sand bath - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

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CERAN® - hot plates

CERAN® an unusual glass ceramic material which is highly resistant to breakage and temperature changes, free from distortion, permeable to ultra-violet light and highly acid-resistant.

Construction: The low mass CERAN® plate is electrically heated over its full surface and is bedded into a stainless steel frame which is mounted on an internally insulated stainless steel housing.

The temperature control system consists of an electronic regulator, activated by a NiCr-Ni temperature sensor (adjustable 50 ... 500°C) installed in a Aluminium pressurecast housing. Warm-up time, e. g. 500°C in 8 minutes. Four feet, which can be individually adjusted to an approx. height of 8 mm to ensure horizontal stability. Type 4A has 4 supporting rod holders (distance 100 mm) In Types SR and EB, the electrical connection between the control housing and the hot plate is protected by a flexible metal hose.

Baureihe A - mit angebautes Regelgehäuse			Series A - with built-in control housing		
Typ Types	Watt	Volt	B x H x T (mm) W x H x D (mm)	CERAN® mm	Gewicht weight
4 A	1500	230 V, 50-60 Hz	440 x 100 x 300	430 x 140	5,0 kg
11 A	2000	230 V, 50-60 Hz	290 x 100 x 410	280 x 280	5,5 kg
22 A	3000	230 V, 50-60 Hz	290 x 100 x 560	280 x 430	6,5 kg
33 A	4400	3x400 V, N+PE, 50-60 Hz	440 x 100 x 560	430 x 430	9,0 kg
44 A	5700	3x400 V, N+PE, 50-60 Hz	590 x 100 x 560	580 x 430	11,5 kg

Baureihe SR - mit separatem Regelgehäuse			Series SR - with separate control housing		
Typ Types	Watt	Volt	B x H x T (mm) W x H x D (mm)	CERAN® mm	Gewicht weight
4 SR	1500	230 V, 50-60 Hz	440 x 100 x 180	430 x 140	5,0 kg
11 SR	2000	230 V, 50-60 Hz	290 x 100 x 290	280 x 280	5,5 kg
22 SR	3000	230 V, 50-60 Hz	290 x 100 x 440	280 x 430	6,5 kg
33 SR	4400	3x400 V, N+PE, 50-60 Hz	440 x 100 x 440	430 x 430	9,0 kg
44 SR	5700	3x400 V, N+PE, 50-60 Hz	590 x 100 x 440	580 x 430	11,5 kg
Die Leitungslänge zwischen Regler und Heizgerät: 1,3 m. Längere Leitungen gegen Mehrpreis The conducting wire between control housing and hot plate measures 1,3 m.					


Baureihe EB - als Einbaugerät für Tischeinbau			Series EB - as built-in units for bench installation			
Typ Types	Watt	Volt	B x H x T (mm) W x H x D (mm)	CERAN® mm	A x B mm	Gewicht weight
11 EB	2000	230 V, 50-60 Hz	290 x 100 x 290	280 x 280	240 x 240	5,5 kg
22 EB	3000	230 V, 50-60 Hz	290 x 100 x 440	280 x 430	390 x 240	6,5 kg
33 EB	4400	3x400 V, N+PE, 50-60 Hz	440 x 100 x 440	430 x 430	390 x 390	9,0 kg
44 EB	5700	3x400 V, N+PE, 50-60 Hz	590 x 100 x 440	580 x 430	390 x 540	11,5 kg
Die Leitungslänge zwischen Regelgehäuse und Heizgerät: 0,65 m. Längere Leitungen gegen Mehrpreis The conducting wire between control housing and hot plate measures 0,65 m.						

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
Further Fabrications

High temperature TITANIUM hot plate up to 600°C - Precision hot plates - Stand hot plates with supporting rod - rapid incinerator - Multi-purpose, extremely deep, heating bath - sand bath - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

The right to make modification serving further technical development, is reserved

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- ✓ Robust
- ✓ Durable
- ✓ also in non-stop operation

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Hot Plates, electric for continuous operation

Electrical hotplates with variable temperature control, with or without a thermostatic regulator. Robust construction for continuous performance. Aluminium alloy hotplate. Because of the asymmetrical, long-term heating system (hotplate corners and edges are subjected to more heat) an even temperature is guaranteed over the entire heating surface. Asbest-free insulation. The housing is constructed from high-grade, stainless steel No. 1.4301, whereby the central section has been additionally sprayed in orange.

Four adjustable feet ensure extra stability. Connecting cable approx. 1,7 m. long, 230 Volt 3200 Watt with Schuko plug.

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Hot plates with wattage power control

Technische Daten Best.-Nr. = Typ	Watt	Volt	Temperatur temperature	A mm	B mm	C mm	D mm	Gewicht weight
HD 0	1800	230	max. 370°C	300	300	312	312	10 kg
HD 1	2200	230	max. 370°C	350	350	358	358	12 kg
HE 1	2400	230	max. 370°C	290	440	308	458	13 kg
HD 2	2850	230	max. 370°C	350	350	364	514	18 kg
HD 3-230	4000	230	max. 370°C	430	580	442	592	25 kg
HD 3-400	4000	3x400, N+PE	max. 370°C	430	580	442	592	25 kg

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Hot plates with wattage power control and thermostaiaic control

Technische Daten Best.-Nr. = Typ	Watt	Volt	Temperatur temperature	A mm	B mm	C mm	D mm	Gewicht weight
HT 01	1000	230	30..110°C	300	300	312	312	11 kg
HT 02	1800	230	50..300°C					
HT 03	1800	230	100..370°C					
HT 11	1150	230	30..110°C	350	350	358	358	13 kg
HT 12	2200	230	50..300°C					
HT 13	2200	230	100..370°C					
TH 11	1650	230	30..110°C	290	440	308	458	14 kg
TH 12	2400	230	50..300°C					
TH 13	2400	230	100..370°C					
HT 21	1800	230	30..110°C	350	500	364	514	19 kg
HT 22	2850	230	50..300°C					
HT 22	2850	230	100..370°C					
HT 31	2000	230	30..110°C	430	580	442	592	26 kg
HT 32-230	4000	230	50..300°C					
HT 32-400	4000	3x400, N+PE	50..300°C					
HT 33-230	4000	230	100..370°C					
HT 33-400	4000	3x400, N+PE	100..370°C	610	160	200	618	11 kg
HB 110	1000	230	30..110°C					
HB 300	2000	230	50..300°C					

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Sand baths, electric for continuous operation

Electrical sandbaths with variable temperature control, with or without a thermostatic regulator. Robust construction for continuous performance. Aluminium alloy hotplate. Because of the asymmetrical, long-term heating system (hotplate corners and edges are subjected to more heat) an even temperature is guaranteed over the entire heating surface. The sandbaths consist of an aluminium alloy hotplate with a screwed-on, sealed stainless steel frame with a usable height of 50 mm. The sand lays direct of the hotplate. The housing is constructed from high-grade, stainless steel No.1.4301, whereby the central section has been additionally sprayed in orange.

4 adjustable feet ensure extra stability. Connecting cable approx. 1,7 m long, 230 Volt/3200 Watt with Schuko plug.

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Sand baths with wattage power control

Technische Daten Best.-Nr. = Typ	Watt	Volt	Temperatur temperature	A mm	B mm	C mm	D mm	Gewicht weight
SD 7	2200	230	max. 370°C	360	360	358	358	13 kg
SD 8	2850	230	max. 370°C	360	360	364	514	21 kg
SD 9-230	4000	230	max. 370°C	440	590	442	592	27 kg
SD 9-400	4000	3x400, N+PE	max. 370°C	440	590	442	592	27 kg

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Sand baths with wattage power control and thermostatic control

Technische Daten Best.-Nr. = Typ	Watt	Volt	Temperatur temperature	A mm	B mm	C mm	D mm	Gewicht weight
ST 71	1150	230	30..110°C	360	360	358	358	14 kg
ST 72	2200	230	50..300°C					
ST 73	2200	230	100..370°C					
ST 81	1800	230	30..110°C	360	500	364	514	21 kg
ST 82	2850	230	50..300°C					
ST 82	2850	230	100..370°C					
ST 91	2000	230	30..110°C	440	590	442	592	28 kg
ST 92-230	4000	230	50..300°C					
ST 92-400	4000	3x400, N+PE	50..300°C					
ST 93-230	4000	230	100..370°C					
ST 93-400	4000	3x400, N+PE	100..370°C					

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Further Fabrications

High temperature TITANIUM hot plate up to 600°C - CERAN hot plates up to 500°C - Electronic DIGITAL hot plates up to 450°C - Precision hot plates - Stand hot plates with supporting rod - Multi-purpose, very deep, heating bath - sand bath - rapid incinerator - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special order

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
Series heating units Hotplate for supporting rod

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Series and supporting rod hot plates

Smoothly finished, cast-iron, hot plate with concave central expansion recess for warp free operation. All housing parts are constructed from stainless steel and each unit has a screw-on fitting for a supporting rod. The appliances have been specifically designed for non-stop use.

Regulation of models SG..., SGK... and SH...

Each unit has a heating output controller that enables a variable adjustment between 10 to 100%.

Regulation of model SGR...

Each unit has an additional, built-in, electronic relay and an extra low-voltage, 3-pole diode output bushing for an optional contact thermometer. When in use, the regulation via the thermometer takes automatic precedence. Without the thermometer, the operation of the SGR... is the same as the SG... model.

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

Further Fabrications

High temperature TITANIUM hot plate up to 600°C - CERAN hot plates up to 500°C - Electronic DIGITAL hot plates up to 450°C - Precision hot plates - Stand hot plates with supporting rod - Multi-purpose, very deep, heating bath - sand bath - rapid incinerator - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

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SVR/E

SVD 95P

- ✓ high-speed heating-up time
- ✓ for single and series incineration of solids and liquid substances

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
RAPID INCINERATOR up to 950°C for single and series incineration of solids and liquid substances

For single and series incineration of solid and liquid substances. The operating temperature of approx. 950°C can be reached within the space of a few minutes. The upper platform has 8 recesses with a diameter of 35 mm which can accommodate porcelain crucibles of up to 50 mm diameter for drying, pre-heating and fuming-off. The uninhibited exposure to oxygen speeds up the incineration process. The inclined position of the crucibles in the lower incinerating bed (which has a length of 400 mm) allows excellent observation and access to the contents during operation. Modern design, easy to operate, reliable, space-saving and of very good value.

Equipment for SV/II (not illustrated):

The upper part and the main housing are made of high-grade stainless steel, the housing comes with an orange finish (RAL 2002). There is a built-in, 120 minute timing switch fitted with a buzzer signal to indicate when the desired time has been reached, or the switch can be placed in the permanently ON position. The appliance is fitted with a safety switch which cuts off the mains supply to the heater immediately the service flap is opened. When the flap is closed, power is restored to the element. The ready-to-use apparatus is delivered complete with a 1.8 m cable and Schuko plug.

Equipment for SVR/E:

As SV/II but with an additional electronic temperature regulation by means of a performance control.

Equipment for SVD 95P:

As SV/II, but with an additional PID temperature regulator up to 999°C, this is controlled by a built-in, NiCr-Ni thermo couple. The pre-selected temperature is shown on the digital display with a constant ACTUAL temperature.

Best-Nr.	SV / II	SVR/E	SVD 95P	Order No.
Spannung	230 - 240 V, 50-60 Hz			voltage
Leistung	2500 Watt			performance
Breite x Tiefe	450 mm x 180 mm			width and depth
Höhe	310 mm			height
max. Tiegel-Ø	50 mm			max. crucible diameter
max. Tiegelhöhe	50 mm			max. crucible height
Gewicht	7,4 kg	8,0 kg	7,8 kg	weight
Best.-Nr. Ersatzheizelement	EH 400	EH 400	EH 95	Order No. spare heating element

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Further Fabrications

High temperature TITANIUM hot plate up to 600°C - CERAN heating plates to 500°C - Electronic DIGITAL hot plates to 450°C - Precision hot plates - Stand hot plates with supporting rod - Multi-purpose, extremely deep, heating bath - sand bath - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

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110°C

240°C

**MULTIPURPOSE
HEATING EQUIPMENT**
Fluid - Oil - Air - Sand Bath
for continuous operation

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W 16 PZ 46-4 PZ 26-7 MB 28/P

110°C 400°C 700°C 240°C

W16, a water bath for temperatures up to 110°C. The container, made of special steel, is compactly built into a container which is also constructed from special steel. The temperature is controlled by a Thermostat which is exact up to $\pm 1,5$ K and more.

A dry-go protection prevents overheating when there is a water shortage. Added protection is also available with the temperature fuse. A deep lying Ring set (Order No. RI 17/E) allows the condensed water to fall back into the bath.

PZ 26-4, which a capacity of 5 litres, can be used as fluid, air and sand bath up to 400°C. As a bath-filling, the combinations of oils, watersoldering tins, sand or other bath-fillings can be used without detrimentally affecting the bath, wh@ is made out of special steel - W.No. 1.4301.

PZ 26-7, designed for temperatures of up to 700°C offers the same possibilities but not, however, as a fluid bath. This apparatus can only be used as a dry air or sand bath. Both baths are also well-suited for use as rotary evaporators.

PZ26-4 and PZ26-7 are both heated with a tubular radiator.

MB 28/P is a fluid bath with an especially large immersion depth of 280 mm, with a bath diameter of 220 mm. The special steel bath container, material No. 1.4301 has no bore holes, a large radius at the bottom and is equipped with a 6-step switchable underfloor heating device switchable performances: 135-220-300-850-1150-1500 Watt. The instrument is equipped with a built-in-thermo-monitor.

PZ26.. and MB28/P are microprocessor controlled by a PID Regulator with a low-mass sensor. The actual value temperature is always on view. To see the temperature setpoint a button has to be pressed. The setpoint temperature is adjustable with an up/down-button.

The housing is constructed from special steel and the base is lacquered. All types of bath are constructed and build for continuous operation, even at maximum temperatures.

Their housings and all component parts will heat up only slightly.

W16, P26-4 u. PZ26-7 in equipped with a screw-on fixture for a support rod of 12 mm dia. MB 28/P for 2 support rod.

Technische Daten Best.-Nr. = Typ		W 16	PZ 26-4	PZ 26-7	MB 28/P	Technical Data Order No. = types
Temperaturvorwahl	°C	30...110	20...400	20...700	20...240	Range of temperature
Anzeige		ohne/without	digital	digital	digital	Display
Leistung	W	1000	1200	1200	1500	Performance
Spannung	V	220...240 V, 50-60 Hz				Voltage
Außenabmessung	mm	220 x 220	310 x 310			Overall dim.
Höhe	mm	160	280	280	435	Height
Gewicht	kg	3,0	6,8	6,8	10,5	Weight

Zubehör	Typ / types	Accessories
für W 16		for W 16
Stativstab 12 mmØ	ST 12	Support rod 12 mm dia.
Ringsatz 175 mmØ	RI 17/E	Ring set 175 mm dia.
für PZ 26..., MB 28/P		for PZ 26..., MB 28/P
Stativstab 12 mmØ	ST 12	Support rod 12 mm dia.
Ringsatz 245 mmØ	RI 24/E	Ring set 245 mm dia.
Reagenzglaseinsatz für 40 Gläser bis 20 mmØ	RG 40	Test glass insert for 40 glasses up to 20 mm dia.
für PZ 26...		for PZ 26...
Lochplatte 245 mmØ für Gefäße mit flachen Böden	LO 90	Perforated plate 245 mm dia. for vessels with flat bottoms
Schale 245 mmØ, 80 mm hoch für Sand o. sonstige Substanzen	SA 24	Cup 245 mm dia., height 80 mm for sand or other substances

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Multiheat
Multi Purpose Heating Equipment
up to 600 °C

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During the construction of the multifunctional heating appliance, all possible laboratory receptacles, regardless of shape or size, were taken into consideration, for heating from the lowest to the highest temperatures. The heating surface of all multifunctional heating appliances consists of a seamless, somewhat conical, high quality stainless steel container with a built-in, insulated tubular heating element in the bottom. The heating conductor is completely insulated and contact-free embedded; it has a spiral form and is soldered to the bottom of the container.

Use as air bath: The low-mass, radiant heater can reach the highest temperature of approx. 800°C in about 5 minutes. The insulated, built-in energy output regulator allows a stage-less calibration of the power from 10...100%, which is equal to a radiated temperature from around 60...600°C.

Use as liquid bath: Because the heating element is completely insulated, the container can be directly filled with e.g. water or oil.


Use as sand bath: with the use of stainless steel, sand bath dishes (which can be delivered) the multifunctional heating appliance can also be used as a sand bath.

Regulation of models STM 50, MU 85, MS 26 and SM...

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Regulation of model MU 85 ER, SMR...

Each unit has an additional, built-in, electronic relay and an extra low-voltage, 3-pole diode output bushing for an optional contact thermometer. When in use, the regulation via the thermometer takes automatic precedence. Without the thermometer, the operation of the units is the same as the MU 85 or SM... model.

												
Techn. Daten Best.-Nr. = Typ	Technical data Ord. No. = type	MU 85	MU85/ER	MS 26/E	STM 50	SM 1 SMR 1	SM 2 SMR 2	SM 3 SMR 3	SM 4 SMR 4	SM 5 SMR 5	SM 6 SMR 6	
Heizstellen	No. of units					1	2	3	4	5	6	
Badinhalt	bath volume	2		5	1	1	2 x 1	3 x 1	4 x 1	5 x 1	6 x 1	L
Breite	Width	225 Ø		310	330	180	355	530	705	880	1055	mm
Tiefe	Depth			310	170	210						
Höhe	Height	280		280	150	225						mm
Badabstand	Distance betw. baths					175						mm
Leistung	power	800		1200	500	500	1000	1500	2000	2500	3000	Watt
Spannung	Voltage	230		230	230	230						Volt
Gewicht	Weight	2,2	2,4	6,9	1,6	2,3 2,4	3,9 4,1	5,5 5,8	7,1 7,5	8,7 9,2	10,3 10,9	kg

Zubehör Accessories	Stativstab Support rod	Ringsatz aus Edelstahl Stainless-steel ring set			Lochplatte aus Edelstahl perforated plate			Sandbadschale sandbath crucible			Halter für Kontaktthermometer Holder for contact thermometer	Kontaktgeberkabel cable f. contact-therm.
Best.-Nr. Order No.	ST 12	RI 14/E	RI 18/E	RI 24/E	LO 14/E	LO 18/E	LO 24/E	SA 14	SA 18	SA 24	HK 3	SK 80/83
SM....	•	•			•			•			nur/only SMR...	nur/only SMR...
STM 50		•			•							
MU 85....	•		•			•			•		nur/only MU 85/ER	nur/only MU 85/ER
MS...	•			•			•			•		

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Further Fabrications

Glass Ceramic Hotplate, High temperature TITANIUM hot plate up to 600°C - Precision hot plates - Stand hot plates with supporting rod - rapid incinerator - Multi-purpose, extremely deep, heating bath - sand bath - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

The right to make modification serving further technical development, is reserved.

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PRÄZITHERM

Precision Hot Plates for continuous operation

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PZ 35

- ✓ Precise temperatures
- ✓ Easy to operate
- ✓ Long-lasting

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Precision Hot Plates for non-stop operation

The heating surfaces of the precision hot plates are composed of an efficient heat conducting Aluminium alloy, and have a smooth anodised finish for easy maintenance. The electrical long-term heating ensures that the temperature remains constant, even on the corners and edges.

The microprocessor controlled digital regulator allows the operating temperature to be preset in 0.1°C stages up to 99.9°C and in 1.0°C stages for temperatures over 100°C. The preset temperature remains programmed, even when the equipment has been switched off. Presetting is only possible by using two fingers; this is to prevent any unintentional alterations that could occur by accidentally touching only one key. The actual temperature is permanently displayed by 14mm high numerals after the keys have been released.

An adjustable safety device against overheating prevents the hotplate and the substance to be heated, from dangerously high temperatures. A removable Plexiglass cover prevents any unauthorised adjustments.

All models have an electronic output control that allows the heating output to be adjusted between 10..100%, this prevents the hotplate exceeding the preset temperature and/or avoids any after-heating. The temperature of all hotplates is measured with a built-in, quick-reaction sensor (PT 100).

All models are equipped with a 3-pole diode output bushing (extra-low voltage), to permit the connection of one of the following appliances:

1. Contact Thermometer – to measure the temperature directly from the heated material and to regulate.
2. Time switch – in order to heat or cool the apparatus at preset intervals.

Model PZ 60 (long shape) is equipped with 5 screw-in sockets for support rods, at intervals of 125 mm.

Special Feature:

Due to the asymmetrical, quick-reaction special heating and thermal sensors as well as output adaptation to the thermal requirements, very precise temperatures with a good uniform heat are obtainable over the whole heating surface.

Technische Daten Best.-Nr. = Typ		PZ 28-1	PZ 28-2	PZ 35	PZ 60	PZ 72		Technical data Order-No. = types
Plattenmaße	mm	200 x 280	200 x 280	350 x 350	610 x 160	430 x 580	mm	Plate size
Temperaturvorwahl	°C	20..110	20..300				°C	Temperature range
Auflösung	K	bis 99,9°C 0,1K, über 100°C 1K					K	Division
Schaltdifferenz	±K	0,1	0,5				±K	Constant temperature
Begrenzer	°C	30..110	50..300				°C	Limiter
Leistung	Watt	500	1100	2200	2000	3300	Watt	Performance
Spannung	Volt	230 V, 50-60 Hz *					Volt	Voltage
Außenabmessung	mm	210 x 300		365 x 380	620 x 215	445 x 610	mm	Overall dimensions
Höhe	mm	135		155		185	mm	Height
Gewicht	kg	7		14	12	22	kg	Weight
* PZ 72 auch für 2x400 V,N+PE lieferbar				* PZ 72 is also available for 2x400 V,N+PE				
Gegen Aufpreis sind alle Heizplatten mit schwarz eloxierter oder mit Teflon beschichteter Heizfläche lieferbar.				For an additional charge, all models can be delivered with an anodised black or Teflon coated heating surface				

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

The microprocessor controlled temperature regulator 2860 SR and 2860 EB are temperature controllers for the precision hotplates. After the desired value has been selected by means of the laminated push buttons, the regulator will hold it constant.

The microprocessor controlled ramp regulator PR 5 SR is a complex time/temperature control for the precision hotplates. Pressing the START key enables the programming of a set of a maximum of 5 ramps (each ramp consists of the ramp-time, constant temperature and the duration of the constant temperature). 3 programme sets can be stored in the non-volatile memory bank.

Further Advantages:


- ✓ A programmed lead time enables the start to be activated up to 99 hours in advance.
- ✓ Easy and comprehensible programming because each temperature and time parameter has its own key. At any time during the running of the programme it is possible to check each programmed value.
- ✓ Output regulator from 0 – 100%
- ✓ Change-over switch for non-stop operation.
- ✓ Programme back-up by EEPROM.
- ✓ Buzzer signal at the end of the programme.

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

Precision Hotplate DIGITAL ELECTRONIC

Temperature Range 20....450°C
with special switching stage for
temperature drop-offs

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- ✓ Continuous operation
- ✓ Precise
- ✓ Robust
- ✓ Reliable

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Precision hot plate to 450°C Type PZ 44

PZ 44 is a hot plate with automatic control for precise temperatures from 20 to 450°C. The required temperature is preset on the digital regulator and remains automatically constant. The actual temperature is continuously displayed with 13 mm high numbers.

The long-lasting electrical heating unit can be used in three heating stages - 825 - 1650 - 3300 Watt. Each stage can be finely tuned from 10 to 100% by an electronic output control to prevent overheating and any after-heating due to excessive power. The heating stages 825 and 3300 Watt, heat the total area of the plate. The special **1650 Watt stage** heats only the right-hand section of the plate. Due to creeping heat, the temperature drops off from the right-hand section to the left-hand plate edge, by approx. 40% (temperature decline).

A low voltage 3-pole socket mounted on the right-hand side of the housing allows the following extras to be connected (however not conditional):

1. Timeswitch, for switching the heating unit on and off at preset times.
2. Contact Thermometer
3. Temperature fuse cut-out

As temperature fuses must have heat contact, a 4 mm diameter hole has been bored on the edge of the right-hand plate for their insertion. The heating unit is automatically cut off from the mains supply after the fuse blows and the contact has reverted to the off position. The fuse has to be renewed once it has blown.

Details and structure of PZ 44

The hot plate is constructed of solid, low-distortion GG 15 casting. Cold housing, also in continuous operation. Material No. 1.4301, from rustproof high grade steel. Centre section is additionally sprayed in orange. Four adjustable feet for extra stability. Sturdy construction, 150 kg max. hot plate load. Ready for connection with 1.5 m. cable.

Best-Nr.	PZ 44/230	PZ 44/400	Order No.
Spannung	230-240 V, 50-60 Hz	2 x 400 V, N+PE	voltage
Leistung	3300 Watt		performance
Plattengröße	440 mm x 290 mm		plate size
Temperaturvorwahl	20°C ... 450°C		Temperature range
Schaltdifferenz	±1 K		Constant temperature
Breite x Tiefe	310 mm x 475 mm		width and depth
Höhe	205 mm		height
Gewicht	26 kg		weight
Zubehör	Best.-Nr. / Order-No.	Accessories	
Kabel mit Stecker für Temperatur-Sicherung, Schaltuhr oder Kontaktthermometer	S 3 P	Cable with plug for temperature-fuse, timeswitch or contact-thermometer	
Temperatur-Sicherungen (± 5°C) 82-128-156-170-182-212-228-254°C (Mindestabnahme 5 Stück)	TS... Wert bitte einsetzen Please ...	Temperature-fuse (± 5°C) 82-128-156-170-182-212-228-254°C	

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Further Fabrications


High temperature TITANIUM hot plate up to 600°C - CERAN heating plates to 500°C - Precision hot plates - Stand hot plates with supporting rod - rapid incinerator - Multi-purpose, extremely deep, heating bath - sand bath - Series heating units - Multipurpose heating equipment - Water- and Oil baths - Development - Special orders.

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PRÄZITHERM
High temperature hot plate of
TITANIUM up to 600 °C
for continuous operation

with PID – Temperature Controller
 or
 Programmer

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of Titanium up to 600°C

PZ 28-3T

This high temperature hot plate reaches its max. temp. of 600°C within about 20 min. The upper plate consists of TITANIUM, overall dim. 280 x 200 mm. The marked area of constant temperature is 230 x 160 mm. Owing to a special insulation the device can be placed even on surfaces sensitive to heat. Metal sheets all around keep the casing temperature to a minimum.

The insulated snap-lid ensures little loss of heat and short times of warming. Between the closed lid (inner surface Ceran®) and the heating plate remain 26 mm, enabling an intensive heating of plain parts at high precise temperatures. The whole casing is made of stainless steel.

The temperature regulator is connected by two plugged cables.

PZ 28-3T is available either with a temperature controller or with a programmer. The temperature controller keeps the device constant to an adjusted temperature within the whole range.

The programmer allows the control of certain programs including phases of heating up, constant temperature and cooling.

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TEMP. CONTROLLER TR 28-3T

This microprocessor controlled temperature regulator is suitable for the high temperature hot plate. After adjustment of the temperature setpoint by the membrane keypad, the PID controller keeps it constant (control deviation 2 K or less).

PROGRAMMER PR 5-3T

The microprocessor controlled ramp regulator offers a comprehensive time/temperature regulation for the high temp. hot plate. Pressing the START key activates a program of up to 5 ramps. Max. 5 ramps (consisting of ramp time, constant temp. and time of constant temp.) can be programmed. 3 program sets can be loaded into the non-volatile memory. A programmed lead time enables the activation of the start up to 99 hours in advance.

Every temperature and time parameter can be selected with its own key. Consequently the setting is easy and clear. It is possible to pass over single ramps by setting the according parameter to '0'. Every set value can be controlled during a program.

Back-up by an EEPROM. The device is equipped with a change-over switch to continuous control; buzzer signal at the program.

Equipment of PR5-3T and TR 28-3T

2-poled main switch; individual adjustment of proportional range, derivative/integrating factor and cycle time; control range limitation (to avoid the setting of values not useful or harmful to an application); contact maker input, actual value correction, junction bush for high temperature hot plate, wattage power control 1...100%, uninterrupted switching.

Both regulators show the actual value constantly. Furthermore the programmer shows the remaining time of a ramp. When the SET key is pressed, the display changes to show the setpoint.

On request the devices are available with a temperature limiter to avoid superheating.