



CORRESPONDING TABLE OF FILTER PAPERS

	WHATMAN	SCHLEICHER & SCHUELL	MACHEREY-NAGEL	SARTORIUS
---	---------	----------------------	----------------	-----------

QUANTITATIVE FILTER PAPERS (ASHLESS)

QT41	N° 41	589/1 black	640w	388
QT42	N° 43	589/2 white	640m	389
QT43	N° 40	589/5 red	640md	392
QT44	N° 44	589/6 green	640dd	390
QT45	N° 42	589/3 blue	640d	391
QT46	-	-	640de	393
QT49	-	589/4 yellow	-	389F
QT51	N° 541	1505	1640w	1388
QT53	N° 540	1506	1640md	1392
QT55	N° 542	1507	1640d	1391

QUALITATIVE FILTER PAPERS

QL01	N° 4	604	617	288 & 1288
QL02	-	597	616	289 & 1289
QL03	N° 2	593	616md	292a & 1292
QL04	N° 6	594	619	290 & 1290
QL05	N° 5	602h	619eh	291 & 1291
QL08	-	-	619de	293

STANDARD FILTER PAPERS

ST60	N° 93	859	-	3hw
ST61	N° 1	595	615	292
ST62	N° 114	503	713	603
ST64	N° 113	520 b II	651	-

GLASS & QUARTZ MICROFIBER FILTERS

FV21	GF/A	GF50	GF-1	GMF 1
FV22	GF/B	GF51	GF-2	GMF 2
FV23	GF/C	GF52	GF-3	GMF 3
FV24	GF/D	GF53	GF-4	GMF 4
FV25	GF/F	GF55	GF-5	GMF 5
FV26	934-AH	GF30	GF-6	GMF 6
FV27	-	GF10	85/90	-
FV29	-	GF6	-	-
FQ30	QM/A	QF20	QF10	QMF