



EURO-SCIENTIFIC FILTER PAPERS



WE GO COOLBLUE



*By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card on the Coolblue site
<https://www.coolblue.be>*

ANALYTICAL QUALITATIVE FILTERS (ash content: 0,06 %)

Ref. QL01 Very fast filtering paper

- Std. Weight : 80 g/m²
- Thickness : 0,21 mm
- Pore size : 15 - 20 µm
- Filtration speed : 10 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL01A0045	4,41	-	-
55	QL01A0055	4,94	-	-
70	QL01A0070	5,54	QL01P0070	27,62
90	QL01A0090	6,63	QL01P0090	27,62
110	QL01A0110	7,76	QL01P0110	29,20
125	QL01A0125	9,96	QL01P0125	30,71
150	QL01A0150	12,15	QL01P0150	32,74
185	QL01A0185	15,42	QL01P0185	36,56
210	QL01A0210	19,88	QL01P0210	42,77

Ref. QL02 Fast filtering paper

- Std. Weight : 88 g/m²
- Thickness : 0,18 mm
- Pore size : 12 - 15 µm
- Filtration speed : 20 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL02A0045	4,41	-	-
55	QL02A0055	4,94	-	-
70	QL02A0070	5,54	QL02P0070	27,62
90	QL02A0090	6,63	QL02P0090	27,62
110	QL02A0110	7,76	QL02P0110	29,20
125	QL02A0125	9,96	QL02P0125	30,71
150	QL02A0150	12,15	QL02P0150	32,74
185	QL02A0185	15,42	QL02P0185	36,56
210	QL02A0210	19,88	QL02P0210	44,07

Ref. QL03 Medium filtering paper

- Std. Weight : 87 g/m²
- Thickness : 0,18 mm
- Pore size : 8 - 12 µm
- Filtration speed : 50 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL03A0045	4,41	-	-
55	QL03A0055	4,94	-	-
70	QL03A0070	5,54	QL03P0070	27,62
90	QL03A0090	6,63	QL03P0090	27,62
110	QL03A0110	7,76	QL03P0110	29,20
125	QL03A0125	9,96	QL03P0125	30,71
150	QL03A0150	12,15	QL03P0150	32,74
185	QL03A0185	15,42	QL03P0185	36,56
210	QL03A0210	19,88	QL03P0210	44,07



WE GO COOLBLUE



**By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card
on the Coolblue site
<https://www.coolblue.be>**

Ref. QL04 Slow filtering paper

- Std. Weight : 80 g/m²
- Thickness : 0,16 mm
- Pore size : 4 - 7 µm
- Filtration speed : 100 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL04A0045	4,41	-	-
55	QL04A0055	4,94	-	-
70	QL04A0070	5,54	QL04P0070	27,62
90	QL04A0090	6,63	QL04P0090	27,62
110	QL04A0110	7,76	QL04P0110	29,20
125	QL04A0125	9,96	QL04P0125	30,71
150	QL04A0150	12,15	QL04P0150	32,74
185	QL04A0185	15,42	QL04P0185	36,56
210	QL04A0210	19,88	QL04P0210	44,07

Ref. QL05 Very slow filtering paper

- Std. Weight : 80 g/m²
- Thickness : 0,16 mm
- Pore size : 2 - 4 µm
- Filtration speed : 200 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL05A0045	4,41	-	-
55	QL05A0055	4,94	-	-
70	QL05A0070	5,54	QL05P0070	27,62
90	QL05A0090	6,63	QL05P0090	27,62
110	QL05A0110	7,76	QL05P0110	29,20
125	QL05A0125	9,96	QL05P0125	30,71
150	QL05A0150	12,15	QL05P0150	32,74
185	QL05A0185	15,42	QL05P0185	36,56
210	QL05A0210	19,88	QL05P0210	44,07



SCIENTIFIC FILTER



WE GO COOLBLUE

By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card on the Coolblue site <https://www.coolblue.be>

ANALYTICAL QUANTITATIVE FILTERS
(ash content: <0,01 %)

Ref. QT41 Very fast filtering paper

- Std. Weight : 85 g/m²
- Thickness : 0,2 mm
- Pore size : 25 - 30 μm
- Filtration speed : 9 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT41A0045	5,51	-	-
55	QT41A0055	6,31	-	-
70	QT41A0070	7,67	QT41P0070	52,35
90	QT41A0090	10,89	QT41P0090	52,35
110	QT41A0110	15,29	QT41P0110	59,54
125	QT41A0125	19,11	QT41P0125	62,91
135	QT41A0135	21,88	-	-
150	QT41A0150	26,24	QT41P0150	75,19
185	QT41A0185	37,18	QT41P0185	102,45
210	QT41A0210	62,81	QT41P0210	135,31
240	QT41A0240	74,85	QT41P0240	163,13

Ref. QT42 Fast filtering paper

- Std. Weight : 85 g/m²
- Thickness : 0,2 mm
- Pore size : 20 - 25 µm
- Filtration speed : 27 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT42A0045	5,51	-	-
55	QT42A0055	6,31	-	-
70	QT42A0070	7,67	QT42P0070	52,36
90	QT42A0090	10,89	QT42P0090	52,36
110	QT42A0110	15,29	QT42P0110	59,54
125	QT42A0125	19,11	QT42P0125	62,91
135	QT42A0135	21,88	-	-
150	QT42A0150	26,24	QT42P0150	75,19
185	QT42A0185	37,18	QT42P0185	102,45
210	QT42A0210	62,81	-	-
240	QT42A0240	74,85	QT42P0240	163,13

Ref. QT43 Medium filtering paper

- Std. Weight : 85 g/m²
- Thickness : 0,2 mm
- Pore size : 14 - 18 µm
- Filtration speed : 55 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT43A0045	5,51	-	-
55	QT43A0055	6,31	-	-
70	QT43A0070	7,67	QT43P0070	52,36
90	QT43A0090	10,89	QT43P0090	52,36
110	QT43A0110	15,29	QT43P0110	59,54
125	QT43A0125	19,11	QT43P0125	62,91
135	QT43A0135	21,88	-	-
150	QT43A0150	26,24	QT43P0150	75,19
185	QT43A0185	37,18	QT43P0185	102,45
210	QT43A0210	62,81	-	-
240	QT43A0240	74,85	QT43P0240	163,13

***By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card
on the Coolblue site
<https://www.coolblue.be>***

Ref. QT44 Slow filtering paper

- Std. Weight : 70 g/m²
- Thickness : 0,16 mm
- Pore size : 7 - 9 µm
- Filtration speed : 140 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT44A0045	5,51	-	-
55	QT44A0055	6,31	-	-
70	QT44A0070	7,67	QT44P0070	52,36
90	QT44A0090	10,89	QT44P0090	52,36
110	QT44A0110	15,29	QT44P0110	59,54
125	QT44A0125	19,11	QT44P0125	62,91
135	QT44A0135	21,88	-	-
150	QT44A0150	26,24	QT44P0150	75,19
185	QT44A0185	37,18	QT44P0185	102,45
210	QT44A0210	62,81	-	-
240	QT44A0240	74,85	QT44P0240	163,13

Ref. QT45 Very slow filtering paper

- Std. Weight : 85 g/m²
- Thickness : 0,17 mm
- Pore size : 2 - 4 µm
- Filtration speed : 195 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT45A0045	5,51	-	-
55	QT45A0055	6,31	-	-
70	QT45A0070	7,67	QT45P0070	52,36
90	QT45A0090	10,89	QT45P0090	52,36
110	QT45A0110	15,29	QT45P0110	59,54
125	QT45A0125	19,11	QT45P0125	62,91
135	QT45A0135	21,88	-	-
150	QT45A0150	26,24	QT45P0150	75,19
185	QT45A0185	37,18	QT45P0185	102,45
210	QT45A0210	62,81	-	-
240	QT45A0240	74,85	QT45P0240	163,13

Other types are available on request.

CELLULOSE EXTRACTION THIMBLES



Reference	Ø ixH (mm)	Packing	Price/unit €	Reference	Ø ixH (mm)	Packing	Price/unit €
ETC16/100	16x100	25	68,83	ETC30/80	30x80	25	67,39
ETC19/90	19x90	25	58,86	ETC30/100	30x100	25	67,39
ETC22/60	22x60	25	74,18	ETC33/80	33x80	25	62,61
ETC22/80	22x80	25	62,61	ETC33/100	33x100	25	68,31
ETC25/80	25x80	25	62,61	ETC33/118	33x118	25	72,12
ETC25/100	25x100	25	67,39	ETC34/130	34x130	25	100,04
ETC26/60	26x60	25	67,39	ETC37/130	37x130	25	81,60
ETC27/80	27x80	25	82,17	ETC41/123	41x123	25	96,77
ETC28/80	28x80	25	67,39	ETC41/150	41x150	25	113,43
ETC28/100	28x100	25	67,39	ETC43/123	43x123	25	120,13

UNIVERSAL INDICATOR PAPER

- Economic range: pH 1-14
- Very quickly and accurately read by a colour scale on the packaging



Reference	Description	Price (€)
PH114R	Universal indicator pH 1-14 (roll of 5 m)	6,50/roll

Reference	Description	Price (€)
8.208.001	Universal indicator pH 1-14 (200 strips)	7,03/200 pcs

Other types and ranges of indicator paper also available at very favourable terms !!!

**Prices in Euro, excl. VAT
Offer gift cards valid until 31/08/2024 !!**



WE GO COOLBLUE

