

EURO-SCIENTIFIC bvba
LERENVELD 12
BE-2547 LINT
TEL. +32 3 489 32 42
FAX +32 3 489 32 73
E-mail info@euro-scientific.be
www.euro-scientific.be

EURO-SCIENTIFIC



Geratherm®

PRODUCT OVERVIEW 2017

Medical Warming Systems



PatientWarmingSystem

UniqueTemp®




UniqueTemp°

The revolutionary PatientWarmingSystem

Geratherm UniqueTemp°

THE CONTROL UNIT | SMALL. CLEVER. INTUITIVE.

article number:	8001140800	
temperature pre-setting:	any connected blanket respectively table pad adjustable separately in steps of 0.1 °C: - OR- and recovery warming blankets for adults: 30 - 40 °C - OR- and recovery warming blankets for children: 30 - 40 °C - OR-operating table pad „Jelly“: 30 - 39 °C	
temperature display:	desired and actual temperature visible on colour display (TFT 3.5")	
operation:	touch-screen menu navigation and selection	
firmware:	updatable to the latest version	
dimensions:	23.2 x 11 x 16.5 cm	
weight:	3.6 kg	
voltage:	input voltage: 100 - 240 V, AC output voltage: 35.5 V, DC	
safety features:	autocheck, temperature measurement between warming system and patient, constant temperature monitoring, autostop in cases of overheating	
connection:	push & pull quick-release coupling (as bus-system)	
classification under directive 93/42 EEC: II b	safety class: I	degree of protection: IPX2



With one Geratherm UniqueTemp° control unit the following elements are usable simultaneously:

4X OR-warming blankets + 1X OR-operating table pad „Jelly“ OR 2X recovery warming blankets

OR-WARMING BLANKETS AND OR-OPERATING TABLE PAD „JELLY“

Hygienic. Cosy. Versatile.



FOR ADULTS

	article	article number	dimensions	weight
	arm/shoulder blanket	8001140202	175 x 32 cm	0.90 kg
	cover for arm/shoulder blanket	8001140203	175 x 32 cm	0.40 kg
	torso blanket	8001140301	52 x 92 cm	0.75 kg
	cover for torso blanket	8001140303	52 x 92 cm	0.20 kg
	leg blanket – right	8001140408	52 x 92 cm	0.75 kg
	cover for leg blanket – right	8001140410	52 x 92 cm	0.20 kg
	leg blanket – left	8001140409	52 x 92 cm	0.75 kg
	cover for leg blanket – left	8001140411	52 x 92 cm	0.20 kg
	operating table pad „Jelly“	8001140600	93 x 49 cm	2.20 kg
	cover for operating table pad „Jelly“	8001140601	93 x 49 cm	0.20 kg

operating voltage: 35.5 V, DC length connecting cable: 70 cm each for ingoing and outgoing leads

FOR CHILDREN

	article	article number	dimensions	weight
	arm/shoulder blanket <i>Paediatric</i>	8001240202	113 x 26 cm	0.60 kg
	cover for arm/shoulder blanket <i>Paediatric</i>	8001240203	113 x 26 cm	0.20 kg
	torso blanket <i>Paediatric</i>	8001240301	71 x 36 cm	0.40 kg
	cover for torso blanket <i>Paediatric</i>	8001240303	71 x 36 cm	0.15 kg
	leg blanket <i>Paediatric</i> – right	8001240408	71 x 36 cm	0.40 kg
	cover for leg blanket <i>Paediatric</i> – right	8001240410	71 x 36 cm	0.15 kg
	leg blanket <i>Paediatric</i> – left	8001240409	71 x 36 cm	0.40 kg
	cover for leg blanket <i>Paediatric</i> – left	8001240411	71 x 36 cm	0.15 kg

UniqueTemp® – state-of-the-art high-performance medical technology for hypothermia prevention.

Its unique microchip technology in combination with its light-weight heat conducting fabric is responsible for the revolutionary characteristics of this elegant electrical patient warming system.

Geratherm UniqueTemp® RECOVERY WARMING BLANKETS

Hygienic. Cosy. Versatile.



FOR ADULTS

article	article number	dimensions	weight
 recovery warming blanket – right	8001150500	120 x 200 cm	2.30 kg
cover for recovery warming blanket – right	8001150502	120 x 200 cm	0.90 kg
 recovery warming blanket – left	8001150501	120 x 200 cm	2.30 kg
cover for recovery warming blanket – left	8001150503	120 x 200 cm	0.90 kg

FOR CHILDREN (AVAILABLE SOON)

article	article number	dimensions	weight
 recovery warming blanket <i>Paediatric</i> – right	8001250500	82 x 112 cm	1.15 kg
cover for recovery warming blanket <i>Paediatric</i> – right	8001250502	82 x 112 cm	0.35 kg
 recovery warming blanket <i>Paediatric</i> – left	8001250501	82 x 112 cm	1.15 kg
cover for recovery warming blanket <i>Paediatric</i> – left	8001250503	82 x 112 cm	0.35 kg

operating voltage: 35.5 V, DC length connecting cable: 130 cm

Geratherm UniqueTemp® ACCESSORIES



article	article number	dimensions	weight
holding claw	8001140900		0.75 kg
extension cable	8001140803	150 cm	0.09 kg
mains power cable	8001120801	250 cm	0.20 kg

GERATHERM MEDICAL AG - Germany

phone: +49 36205 98-0 • fax: +49 36205 98-116 • info@geratherm.com • www.geratherm.com