



Dry warming of baby
bottles and syringes

**5 years
warranty**

Barkey clinitherm baby

Babies are the most sensitive patients. Hygiene is very deciding when it comes to warming of baby nutrition. Barkey clinitherm baby warms baby bottles safe, dry and without the risk of contamination with germs.

Benefits at a glance

- Intended for the safe warming of nutrition for babies
- High hygienic standards due to dry warming
- Fits for different bottles sizes and diameters
- Adapters for the use of different bottles/ syringes and diameters (optional)
- Available in different sizes
- Infinite adjustment of temperature in between the range of +20°C to +60°C
- Safety due to rounded edges
- Easy cleaning through wipe disinfection
- Dishwasher proof adapters
- Electronic and mechanical overtemperature alarm

Barkey clinitherm baby

Safety

Electronic and mechanic overtemperature alarms almost eliminate the overheating of baby nutrition. Thereby valuable ingredients of the milk survive.

Furthermore the device is of low noise level which enables the Barkey clinitherm baby to be used in mother-baby rooms.



Product	Barkey clinitherm baby
Capacity	from 1 to 25 bottles
Form	round
Immersion depth	up to 120 mm
Reducing inserts	in accordance to customer's needs
Dimensions	according to model
Casing	aluminium, baked varnish
Warmer	aluminium, oxidized
Weight	according to model
Temperature	ambient temperature to +60°C

Overtemperature protection	electronic +65°C, mechanical +73°C
Mains supply	230 V 50/60 Hz
Power consumption	according to model
Current consumption	according to model
Current fuse	according to model
Protection class	I, protective earthing
Storage/Transport	-20°C to +50°C
Certification	CE
Barkey is certified	DIN EN ISO 13485

* The basic model can be altered to suit the customer's individual requirements.