



## EURO-SCIENTIFIC FILTER PAPERS



**WE GO COOLBLUE**



*By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card on the Coolblue site  
<https://www.coolblue.be>*

### ANALYTICAL QUALITATIVE FILTERS (ash content: 0,06 %)

Ref. QL01 Very fast filtering paper

- Std. Weight : 80 g/m<sup>2</sup>
- Thickness : 0,21 mm
- Pore size : 15 - 20 µm
- Filtration speed : 10 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL01A0045	3,83	-	-
55	QL01A0055	4,28	-	-
70	QL01A0070	4,80	QL01P0070	-
90	QL01A0090	5,78	QL01P0090	23,50
110	QL01A0110	6,72	QL01P0110	24,83
125	QL01A0125	8,65	QL01P0125	26,17
150	QL01A0150	10,52	QL01P0150	27,85
185	QL01A0185	13,40	QL01P0185	31,10
210	QL01A0210	17,27	QL01P0210	-

**Ref. QL02    Fast filtering paper**

- Std. Weight           : 88 g/m<sup>2</sup>
- Thickness             : 0,18 mm
- Pore size             : 12 - 15 µm
- Filtration speed     : 20 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL02A0045	3,83	-	-
55	QL02A0055	4,28	-	-
70	QL02A0070	4,80	QL02P0070	23,50
90	QL02A0090	5,78	QL02P0090	23,50
110	QL02A0110	6,72	QL02P0110	24,83
125	QL02A0125	8,65	QL02P0125	26,17
150	QL02A0150	10,52	QL02P0150	27,85
185	QL02A0185	13,40	QL02P0185	31,13
210	QL02A0210	17,27	QL02P0210	37,50

**Ref. QL03    Medium filtering paper**

- Std. Weight           : 87 g/m<sup>2</sup>
- Thickness             : 0,18 mm
- Pore size             : 8 - 12 µm
- Filtration speed     : 50 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL03A0045	3,83	-	-
55	QL03A0055	4,28	-	-
70	QL03A0070	4,80	QL03P0070	23,50
90	QL03A0090	5,78	QL03P0090	23,50
110	QL03A0110	6,72	QL03P0110	24,83
125	QL03A0125	8,65	QL03P0125	26,17
150	QL03A0150	10,52	QL03P0150	27,85
185	QL03A0185	13,40	QL03P0185	31,10
210	QL03A0210	17,27	QL03P0210	37,50



**WE GO COOLBLUE**



**By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card  
on the Coolblue site  
<https://www.coolblue.be>**

**Ref. QL04**    Slow filtering paper

- Std. Weight           : 80 g/m<sup>2</sup>
- Thickness             : 0,16 mm
- Pore size             : 4 - 7 µm
- Filtration speed     : 100 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL04A0045	3,83	-	-
55	QL04A0055	4,28	-	-
70	QL04A0070	4,80	QL04P0070	23,50
90	QL04A0090	5,78	QL04P0090	23,50
110	QL04A0110	6,72	QL04P0110	24,83
125	QL04A0125	8,65	QL04P0125	26,17
150	QL04A0150	10,52	QL04P0150	27,85
185	QL04A0185	13,40	QL04P0185	31,10
210	QL04A0210	17,27	QL04P0210	37,50

**Ref. QL05**    Very slow filtering paper

- Std. Weight           : 80 g/m<sup>2</sup>
- Thickness             : 0,16 mm
- Pore size             : 2 - 4 µm
- Filtration speed     : 200 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QL05A0045	3,83	-	-
55	QL05A0055	4,28	-	-
70	QL05A0070	4,80	QL05P0070	23,50
90	QL05A0090	5,78	QL05P0090	23,50
110	QL05A0110	6,72	QL05P0110	24,83
125	QL05A0125	8,65	QL05P0125	26,17
150	QL05A0150	10,52	QL05P0150	27,85
185	QL05A0185	13,40	QL05P0185	31,10
210	QL05A0210	17,27	QL05P0210	37,50



SCIENTIFIC FILTER



**WE GO COOLBLUE**

*By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card on the Coolblue site <https://www.coolblue.be>*

**ANALYTICAL QUANTITATIVE FILTERS**  
(ash content: <0,01 %)

Ref. QT41 Very fast filtering paper

- Std. Weight : 85 g/m<sup>2</sup>
- Thickness : 0,2 mm
- Pore size : 25 - 30 μm
- Filtration speed : 9 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT41A0045	4,78	-	-
55	QT41A0055	5,47	-	-
70	QT41A0070	6,65	QT41P0070	44,55
90	QT41A0090	9,45	QT41P0090	44,55
110	QT41A0110	13,25	QT41P0110	50,67
125	QT41A0125	16,58	QT41P0125	53,53
135	QT41A0135	18,98	-	-
150	QT41A0150	22,77	QT41P0150	63,95
185	QT41A0185	32,25	QT41P0185	87,18
210	QT41A0210	54,53	-	-
240	QT41A0240	64,93	QT41P0240	138,82

**Ref. QT42    Fast filtering paper**

- Std. Weight            : 85 g/m<sup>2</sup>
- Thickness             : 0,2 mm
- Pore size              : 20 - 25 µm
- Filtration speed     : 27 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT42A0045	4,78	-	-
55	QT42A0055	5,47	-	-
70	QT42A0070	6,65	QT42P0070	44,55
90	QT42A0090	9,45	QT42P0090	44,55
110	QT42A0110	13,25	QT42P0110	50,67
125	QT42A0125	16,58	QT42P0125	53,53
135	QT42A0135	18,98	-	-
150	QT42A0150	22,77	QT42P0150	63,95
185	QT42A0185	32,25	QT42P0185	87,18
210	QT42A0210	54,53	-	-
240	QT42A0240	64,93	QT42P0240	138,82

**Ref. QT43    Medium filtering paper**

- Std. Weight            : 85 g/m<sup>2</sup>
- Thickness             : 0,2 mm
- Pore size              : 14 - 18 µm
- Filtration speed     : 55 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT43A0045	4,78	-	-
55	QT43A0055	5,47	-	-
70	QT43A0070	6,65	QT43P0070	44,58
90	QT43A0090	9,45	QT43P0090	44,58
110	QT43A0110	13,25	QT43P0110	50,67
125	QT43A0125	16,58	QT43P0125	53,53
135	QT43A0135	18,98	-	-
150	QT43A0150	22,77	QT43P0150	63,95
185	QT43A0185	32,25	QT43P0185	87,18
210	QT43A0210	54,53	-	-
240	QT43A0240	64,93	QT43P0240	138,82

***By purchase of each € 60 net Euro-Scientific offers you a € 10 gift card  
on the Coolblue site  
<https://www.coolblue.be>***

Ref. QT44    Slow filtering paper

- Std. Weight           : 70 g/m<sup>2</sup>
- Thickness             : 0,16 mm
- Pore size             : 7 - 9 µm
- Filtration speed     : 140 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT44A0045	4,78	-	-
55	QT44A0055	5,47	-	-
70	QT44A0070	6,65	QT44P0070	44,55
90	QT44A0090	9,45	QT44P0090	44,58
110	QT44A0110	13,25	QT44P0110	50,67
125	QT44A0125	16,58	QT44P0125	53,53
135	QT44A0135	18,98	-	-
150	QT44A0150	22,77	QT44P0150	63,95
185	QT44A0185	32,25	QT44P0185	87,18
210	QT44A0210	54,53	-	-
240	QT44A0240	64,93	QT44P0240	138,82

Ref. QT45    Very slow filtering paper

- Std. Weight           : 85 g/m<sup>2</sup>
- Thickness             : 0,17 mm
- Pore size             : 2 - 4 µm
- Filtration speed     : 195 sec

Ø mm	Plain disc	€/100 pieces	Folded disc	€/100 pieces
45	QT45A0045	4,78	-	-
55	QT45A0055	5,47	-	-
70	QT45A0070	6,65	QT45P0070	44,55
90	QT45A0090	9,45	QT45P0090	44,55
110	QT45A0110	13,25	QT45P0110	50,67
125	QT45A0125	16,58	QT45P0125	53,53
135	QT45A0135	18,98	-	-
150	QT45A0150	22,77	QT45P0150	63,95
185	QT45A0185	32,25	QT45P0185	87,18
210	QT45A0210	54,53	-	-
240	QT45A0240	64,93	QT45P0240	138,82

***Other types are available on request.***

## CELLULOSE EXTRACTION THIMBLES




Reference	Ø ixH (mm)	Packing	Price/unit €	Reference	Ø ixH (mm)	Packing	Price/unit €
ETC16/100	16x100	25	59,72	ETC30/80	30x80	25	58,47
ETC19/90	19x90	25	51,08	ETC30/100	30x100	25	58,47
ETC22/60	22x60	25	64,40	ETC33/80	33x80	25	54,33
ETC22/80	22x80	25	54,33	ETC33/100	33x100	25	59,30
ETC25/80	25x80	25	54,33	ETC33/118	33x118	25	62,58
ETC25/100	25x100	25	58,47	ETC34/130	34x130	25	86,80
ETC26/60	26x60	25	58,47	ETC37/130	37x130	25	70,82
ETC27/80	27x80	25	71,30	ETC41/123	41x123	25	83,97
ETC28/80	28x80	25	58,47	ETC41/150	41x150	25	98,38
ETC28/100	28x100	25	58,47	ETC43/123	43x123	25	104,22

## UNIVERSAL INDICATOR PAPER

- Economic range: pH 1-14
- Very quickly and accurately read by a colour scale on the packaging



Reference	Description	Price (€)
8.208.010	Universal indicator pH 1-14 (roll of 5 m)	5,57/roll



Reference	Description	Price (€)
8.208.001	Universal indicator pH 1-14 (200 strips)	6,05/200 pcs

*Other types and ranges of indicator paper also available at very favourable terms !!!*

**Prices in Euro, excl. VAT  
Offer gift cards valid until 31/08/2021 !!**



**WE GO COOLBLUE**

