

MAILING

"GRAM" REFRIGERATORS AND FREEZERS

BIOBASIC SERIE:

The BIOBASIC range is available as a refrigerator or freezer as underbench (125L), intermediate (218L) or full height (346L) and (600L) models, all equipped with wire shelves as standard. With temperature alarms, internal and external ATEX, no cold walls and a user-friendly interface, the BIOBASIC range offers an exceptional selection of features in a very small package. All refrigerators 210/310/410 are also available with glass doors.



REFRIGERATORS:

Technical specifications	RR210	RR310	RR410	RR600
temperature range	+2/+15°C	+2/+15°C	+2/+15°C	+2/+15°C
ext. dimensions in mm	595x640x831	595x640x1220	595x640x1876	700x895x2125
net volume in litres	104	189	312	536
material interior/exterior	PS/white lacquered steel	PS/white lacquered steel	PS/white lacquered steel	Stainless steel/white lacquered steel
refrigerant	R600a	R600a	R600a	R290
ATEX-marking	II 3G Ex ec nC ic IIB T6 Gc DTI 22ATEX0252X			
art. nr. closed door	782530001	783530001	784530001	786000001
price excl. VAT in €	1630,-	1800,-	1960,-	2930,-
art. nr. window door	782540001	783540001	784540001	-
price excl. VAT in €	1895,-	2065,-	2225,-	-

FREEZERS :

Technical specifications	RF210	RF310	RF410	RF600
temperature range	-25/-5°C	-25/-5°C	-25/-5°C	-25/-5°C
ext. dimensions in mm	595x640x831	595x640x1220	595x640x1876	700x895x2125
net volume in litres	104	189	312	536
material interior/exterior	PS/white lacquered steel	PS/white lacquered steel	PS/white lacquered steel	Stainless steel/white lacquered steel
refrigerant	R600a	R600a	R600a	R290
ATEX-marking	II 3G Ex ec nC ic IIB T3 Gc/DTI 22ATEX0252X			
art. Nr.	782550001	783550001	784550001	786020001
price excl. VAT in €	1740,-	1910,-	2070,-	3130,-

Additional options and information available upon simple request!

BIOCOMPACT II SERIE:

These compact refrigerators and freezers are suitable for a wider range of biostorage purposes where high reliability is a must. The BioCompact II range provides you with significantly better performance than any other cabinets in this segment, on account of a refrigeration system that does away with cold walls that can damage the contents and a unique air distribution system for storing ordinary biomaterial under consistent, stable conditions. The design of the BioCompact II range ensures you lots of ways to deal with individual storage needs, thanks to the versatile interior layout and a comprehensive selection of fixtures and fittings. The small footprint of this design also makes it the perfect biostorage unit for use in confined spaces.



REFRIGERATORS:

Technical specifications	RR210	RR310	RR410	RR610
temperature range	+2/+20°C	+2/+20°C	+2/+20°C	+2/+20°C
ext. dimensions in mm	595x640x801	595x640x1190	595x640x1776	695x875x1874
net volume in litres	104	189	312	419
material interior/exterior	PS/white or stainless steel	PS/white or stainless steel	PS/white or stainless steel	PS/white or stainless steel
refrigerant	R600a	R600a	R600a	R600a
ATEX-marking	II 3G Ex ec nC ic IIB T6 Gc DTI 22ATEX0251X			
art. nr. white	IIRR210LG	IIRR310LG	IIRR410LG	IIRR610LG
price excl. VAT in € white	2160,-	2325,-	2490,-	3310,-
art. nr. stainless steel	IIRR210CG	IIRR310CG	IIRR410CG	IIRR610CG
price excl. VAT in € stainless steel	2270,-	2435,-	2600,-	3410,-

FREEZERS :

Technical specifications	RF210	RF310	RF410	RF610
temperature range	-25/-5°C	-25/-5°C	-25/-5°C	-25/-5°C
ext. dimensions in mm	595x640x801	595x640x1190	595x640x1776	695x875x1874
net volume in litres	104	189	312	419
material interior/exterior	PS/white or stainless steel	PS/white or stainless steel	PS/white or stainless steel	PS/white or stainless steel
refrigerant	R600a	R600a	R600a	R290
ATEX-marking	II 3G Ex ec nC ic IIB T3 Gc/DTI 22ATEX0252X			
art. nr. white	IIRF210LG	IIRF310LG	IIRF410LG	IIRF610LG
price excl. VAT in € white	2275,-	2440,-	2605,-	3615,-
art. Nr. stainless steel	IIRF210CG	IIRF310CG	IIRF410CG	IIRF610CG
price excl. VAT € stainless steel	2385,-	2550,-	2715,-	3715,-

Additional options and information available upon simple request!

BIO ULTRA SERIE:

The BioUltra UL570 is an ultra-low freezer designed to safely store materials down to -86°C. By using very high-quality components, we ensure that the BioUltra are able to maintain safe and stable interior conditions when operating at -86°C. The cabinets from this line are available with a standard air-cooled condenser or optionally with a water- and air-cooled hybrid condenser. Customisable layout allows users to adapt the cabinet to current storage requirement, 5 stainless steel perforated shelves are included as standard.



Technical specifications	UL570	
temperature range	-86/-60 °C	
ext. dimensions in mm	856x979x1986	
brut volume in litres	570	
material interior/exterior	white lacquered steel/ white lacquered steel	white lacquered steel /stainless steel
refrigerant	Stage 1: R290 / Stage 2: R170/601	
ATEX-marking	II 3G Ex ec nC ic IIB T4 Gc/DTI 22ATEX0248X	
price excl. VAT in €	11450,-	11900,-

EX GUARD SERIE:

The ExGuard is designed for safe temperature controlled storage of items potentially emitting noxious, malodorous and/or explosive atmospheres. At the core, the ExGuard is envisioned to facilitate a safer work environment for users.

ExGuard increases work environment safety and allows placement in work environments by letting your existing ventilation system extract the hazardous atmospheres from the cabinet interior prior to you being able to access the interior.

As with all other BioLine refrigerators, the ExGuard has internal and external ATEX Category 3 Zone 2 compliance, thereby keeping the work environment safety as high as possible.



Technical specifications	ER600W	
temperature range	-2/20°C	
ext. dimensions in mm	820x789x1996	
net volume in litres	486	
materiaal intern/extern	stainless steel/white lacquered steel	stainless steel/stainless steel
refrigerant	R290	
ATEX-marking	II 3G Ex ec nC ic IIB T1 Gc/ DTI 22ATEX0249X	
price excl. VAT in €	8070,-	8320,-

Additional options and information available upon simple request!

SHIPPING COSTS WILL BE CALCULATED SEPARATELY!!