

Euro-Scientific News

LABINCO LABORATORY INSTRUMENTS

■ BASIC MAGNETIC STIRRERS

Magnetic stirrers with heating

- analogue
- aluminium heating plate
- delivered with magnetic stirring bar
- CE-certified
- 220 - 230V/50 - 60Hz
- warranty: 24 months



	L81	L32
Dimensions of the heating plate (mm)	Ø 145	Ø 195
Stirring capacity (litres water)	10	20
Temperature range (°C)	50 - 325	50 - 325
Power (W)	600	600
Stirring speed (rpm)	50 - 1250	50 - 1250
Dimensions (mm)	172 x 286 x (H)125	200 x 250 x (H)135
REFERENCE	81000	32000
PRICE (€)	346,50	367,95
PROMO (€)	322,25	342,19

Digital versions available on request!

MULTI-POINT MAGNETIC STIRRERS

Multi-point magnetic stirrers with heating

- temperature range: 25-350 °C (digital)
- speed range: 50-1250 rpm (analogue)
- aluminium heating plates Ø 120 mm
- 4 positions
- up to 5 litres per stirring point
- delivered with 4 magnetic stirring bars
- 220 - 230V/50 - 60Hz
- CE-certified
- warranty: 24 months



	1 x 4 stirring points	2 x 2 stirring points
Type	LD814	LD844
Dimensions (mm)	550 x 245 x (H)140	345 x 385 x (H)185
REFERENCE	81400	84400
PRICE (€)	1720,95	1720,95
PROMO (€)	1600,48	1600,48

1 x 6 or 2 x 3 stirring points available on request!
Option: plates made of stainless steel or PTFE-coated

VORTEX MIXER

Vortex L46

- strong vibrating action, for test tubes Ø 7 mm up to 35 mm
- "Press-to-Mix" and continuous action
- speed range: 50-1250 rpm
- amplitude: 4 mm
- dimensions: 135x120x(H)135 mm
- CE-certified
- warranty: 24 months


REFERENCE	46000
PRICE (€)	280,50
PROMO (€)	260,87



■ ROLLING BENCHES

Rolling benches with strong motor to keep samples in motion

- one roll is driven
- delivered with 2 rolls, 1 fixed and 1 adjustable
- roll Ø 46 mm
- for up to 5 rolls
- with digital timer up to 99h 59m 59s
- analogue speed controller: 1-240 rpm
- 230 V
- CE-certified
- warranty: 24 months



REFERENCE	ROLL LENGTH (mm)	PRICE (€)	PROMO (€)	REF. EXTRA ROLL	PRICE (€)	PROMO (€)
20935	350	1791,90	1666,47	35209	84,15	78,26
20945	450	1826,55	1698,69	45209	92,40	85,93
20955	550	1882,65	1750,86	55209	100,65	93,60
20965	650	1947,00	1810,71	65209	108,90	101,28
20975	750	2009,70	1869,02	75209	117,15	108,95

Rolling benches with timer, memory function and re-start function available on request!

■ TEST-TUBE ROTATOR

Test-tube rotator L28

- test-tube rotator for serial rotation of test-tubes 4-40 mm
- speed 10 rpm, not adjustable
- angle of rotation can be set between 27° and 90°
- dimensions: 200x200x(H)220 mm
- CE-certified
- warranty: 24 months

REFERENCE	28000
PRICE (€)	471,90
PROMO (€)	438,87

Disks for rotator L28

Disk with 40 clamps for tubes Ø 4-6 mm

REFERENCE	28006
PRICE (€)	107,25
PROMO (€)	99,74

Disk with 30 clamps for tubes Ø 6-9 mm

REFERENCE	28009
PRICE (€)	105,60
PROMO (€)	98,21


Disk with 20 clamps for tubes Ø 16-19 mm

REFERENCE	28030
PRICE (€)	90,75
PROMO (€)	84,40

Drum with 60 holes Ø 17,5 mm

REFERENCE	28085
PRICE (€)	155,10
PROMO (€)	144,24

Other disks available on request!



**Prices, excl. VAT and recupel,
valid until 30/05/2024**