



Thawing and warming of
plasma and blood products
– safe and reliable

Barkey plasmatherm

The Barkey plasmatherm thaws FFP bags and warms blood, erythrocyte concentrates (EC), cryopreserved preparations, cryopreserved stem cells (HPC) and infusion solutions.

Benefits at a glance

- Thaws 4 bags at a time (in case of emergency up to 8)
- Simple and intuitive application
- High standard of hygiene due to dry warming (no direct contact to water)
- Smooth surfaces and rounded edges
- Easy to clean through wipe disinfection
- Electronic documentation (optional)
- 2 leak detection sensors = no risk of contamination or infection
- Shortcuts for plasma, blood, HPC and user program
- Large display with plain text
- Change of water just once a year
- Cushions conform to bags



Barkey plasmatherm

Accessories (optional)

- Electronic documentation software (Ethernet connection only)
- Barcode reader
- Log printer

Hygiene

- No risk of contamination due to dry thawing and warming (no direct contact to water)
- Easy to clean through wipe disinfection
- Corrosion-resistant

Safety

- Electronic and mechanical overtemperature alarm
- Adjustable temperature range from +37°C to +45°C
- Visual and acoustic signal after heating duration is complete or in case of error
- Visual inspection due to clear cushions and lid
- 2 leak detection sensors



Product	Barkey plasmatherm
Power supply	220-240 V 50/60 Hz or 100-120 V 50/60 Hz
Set temperature	+37°C to max. +45°C; US: +37°C
Dimensions	340 x 600 x 320 mm (WxDxH)
Weight without filling	18 kg
Filling quantity	9 l

Interfaces/Devices (optional)	Barcode scanner, Log printer, Ethernet port
MDD classification	Ila
Safety class	I, protective earthing
Certification	CE 0123
Barkey is certified	DIN EN ISO 13485